

CARTPT Polyclonal Antibody

Catalog number: 13212-1-AP

Size: 28 µg/150 µl

Source: Rabbit

Isotype: IgG

Synonyms:

CARTPT; CART, CART

prepropeptide, CARTPT



SH-SY5Y cells were subjected to SDS PAGE followed by western blot with 13212-1-AP(CARTPT antibody) at dilution of 1:800

Applications

Tested applications:	ELISA, WB
Species specificity:	Human,Mouse,Rat; other species not tested.
Calculated CARTPT MW:	116aa,13 kDa
Observed CARTPT MW:	4-14kd
Positive WB detected in	SH-SY5Y cells, human adrenal gland tissue
Recommended dilution:	WB: 1:200-1:2000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation

Immunogen information

Immunogen:	Ag4004
GenBank accession number:	BC029882
Gene ID (NCBI):	9607
Full name:	CART prepropeptide

Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3. Store at -20°C.