

RRP9 Polyclonal Antibody

Catalog number: 10311-1-AP

Size: 20 µg/150 µl

Source: Rabbit

Isotype: IgG

Synonyms:

RRP9; RNU3IP2; U3-55K

Applications

Tested applications: ELISA
Species specificity: Human, mouse, rat; other species not tested.
Calculated RRP9 MW: 52kd
Observed RRP9 MW:

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation, FC = Flow cytometry

Immunogen information

Immunogen: Ag0319
GenBank accession number: BC001113
Gene ID (NCBI): 9136
Full name: Ribosomal RNA processing 9, small subunit (SSU) processome component, homolog (yeast)

Product information

Purification method: Antigen affinity purification
Storage: PBS with 0.1% sodium azide and 50% glycerol pH 7.3. Store at -20°C.