

GDAP1 Polyclonal Antibody

Catalog number: 13152-1-AP


Size: 66 µg/150 µl

Source: Rabbit

Isotype: IgG

Synonyms:

GDAP1; CMTRIA, GDAP1



human brain tissue were subjected to SDS PAGE followed by western blot with 13152-1-AP(GDAP1 antibody) at dilution of 1:500

Applications

Tested applications:	ELISA, WB
Species specificity:	Human, Mouse; other species not tested.
Calculated GDAP1 MW:	358aa,41 kDa
Observed GDAP1 MW:	36-41 kDa
Positive WB detected in	Human brain tissue, fetal human brain tissue, mouse brain tissue
Recommended dilution:	WB: 1:200-1:2000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation

Immunogen information

Immunogen:	Ag3769
GenBank accession number:	BC024939
Gene ID (NCBI):	54332
Full name:	Ganglioside-induced differentiation-associated protein 1

Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.1% sodium azide and 50% glycerol pH 7.3. Store at -20°C.