

## BLNK Polyclonal Antibody

Catalog number: 10855-1-AP

Size: 20 µg/150 µl


Source: Rabbit

Isotype: IgG


Synonyms:

BLNK; BASH; BLNK-S; LY57;

MGC111051; SLP-65; SLP65



Raji cells were subjected to SDS PAGE followed by western blot with 10855-1-AP (BLNK Antibody) at dilution of 1:300



Immunohistochemistry of paraffin-embedded human pancreas cancer using 10855-1-AP (BLNK Antibody) at Dilution 1:100 (under 10x lens)

### Applications

Tested applications:	ELISA, WB, IHC
Species specificity:	Human, mouse, rat; other species not tested.
Positive WB detected in:	Raji cells
Calculated BLNK MW:	50kd
Observed BLNK MW:	68kd, 70kd
Positive IHC detected in:	Human pancreas cancer
Recommended dilution:	WB: 1:200-1:1000 IHC: 1:50-1:200

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation, FC = Flow cytometry

### Immunogen information

Immunogen:	Ag1290
GenBank accession number:	BC018906
Gene ID (NCBI):	29760
Full name:	B-cell linker

### Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.1% sodium azide and 50% glycerol pH 7.3. Store at -20°C.