



Cal-Pak Reference Chart



800.345.2100 • 800.345.5232 fax
PO Box 99 • Swedesboro, NJ 08085
www.thomassci.com • value@thomassci.com



Vendor Name	Thomas part #	Category	Description	Manu Part #	Cal-Pak* #	Thomas	Recommended Weight #1			Recommended Weight #2			Recommended Weight #3		
							Size	Part #	Thomas	Size	Part #	Thomas	Size	Part #	Thomas
A&D ENGINEERING	1226U35	BALANCES	BALANCE COMPACT 400g x 0.01g 120VAC	EJ-410			200g	7016-1W	1207W25	20g	7020-1W	0867B41			
A&D ENGINEERING	1226U41	BALANCES	BALANCE COMPACT 6000g x 0.1g 120VAC	EJ-6100			5kg	7010-1W	1207W21	200g	7016-1W	1207W25			
A&D ENGINEERING	1326B51	BALANCES	BALANCE 1200 X .1G	EK-1200i			1kg	7013-1W	1207W23	200g	7016-1W	1207W25			
A&D ENGINEERING	1326B43	BALANCES	BALANCE 120 X .01 G	EK-120i			100g	7017-1W	0867B39	20g	7020-1W	0867B41			
A&D ENGINEERING	1326K45	BALANCES	BALANCE 200 X .01G	EK-200i			200g	7016-1W	1207W25	20g	7020-1W	0867B41			
A&D ENGINEERING	1326B55	BALANCES	BALANCE 3000 X .1G	EK-3000i			2kg	7012-1W	1207W22	200g	7016-1W	1207W25			
A&D ENGINEERING	1326B47	BALANCES	BALANCE 300 X .01G	EK-300i			200g	7016-1W	1207W25	20g	7020-1W	0867B41			
A&D ENGINEERING	1326B49	BALANCES	BALANCE 600 X .01G	EK-600i			500g	7014-1W	1207W24	20g	7020-1W	0867B41			
A&D ENGINEERING	1226T78	BALANCES	BALANCE COMPACT 6000g x 0.1g 110V	EK-6100i			5kg	7010-1W	1207W21	200g	7016-1W	1207W25			
A&D ENGINEERING	1226U16	BALANCES	SCALE BENCH 60 x 0.02lb NTEP 110VAC	FG-30KBMN			50g	7018-1W	0867B40	5g	7022-1W	0867B43			
A&D ENGINEERING	1226T63	BALANCES	BALANCE TOPLOAD 1220g x 0.01g 110V	FX-1200i	7263-1W	1228U04									
A&D ENGINEERING	1226T69	BALANCES	BALANCE IP65 1220g x 0.01g 110V	FX-1200iWP	7263-1W	1228U04									
A&D ENGINEERING	1226T70	BALANCES	BALANCE IP65 2200g x 0.01g 110V	FX-2000iWP	7265-1W	1228U06									
A&D ENGINEERING	1226T68	BALANCES	BALANCE IP65 320g x 0.001g 110V	FX-300iWp	7258-1W	1228T99									
A&D ENGINEERING	1326C53	BALANCES	Toploader Balance 310gx0.001g w/Breeze Break	GF-300	7258-1W	1228T99									
A&D ENGINEERING	1326C62	BALANCES	BALANCE 3100 X 0.01G	GF-3000	7265-1W	1228U06									
A&D ENGINEERING	1326C65	BALANCES	BALANCE 6100 X 0.1G	GF-6000		5kg	7010-1W	1207W21	200g	7016-1W	1207W25				
A&D ENGINEERING	1226T44	BALANCES	BALANCE ANALYTICAL 120g x 0.1mg 120V	GH-120	7255-1W	1228T96									
A&D ENGINEERING	1226T48	BALANCES	BALANCE ANALYTIC 250/101g x 0.01/0.1mg 120V	GH-252		200g	7016-1W	1207W25							
A&D ENGINEERING	1327E01	BALANCES	BALANCE 21000X.1G GP20K	GP-20K		20kg	7008-1W	1207W19	2kg	7012-1W	1207W22	200g	7016-1W	1207W25	
A&D ENGINEERING	1327E03	BALANCES	BALANCE 31000X.1G GP30K	GP-30K		20kg	7008-1W	1207W19	2kg	7012-1W	1207W22	200g	7016-1W	1207W25	
A&D ENGINEERING	1327E05	BALANCES	BALANCE 41000X.5G GP40K	GP-40K		20kg	7008-1W	1207W19	1kg	7013-1W	1207W23				
A&D ENGINEERING	1325B03	BALANCES	Analytical Balance 210gx0.1mg w/RS-232C	GR-200	7257-1W	1228T98									
A&D ENGINEERING	1325B01	BALANCES	Analytical Balance 42g/210gx0.1mg/0.1mg w/ RS-232C	GR-202		200g	7016-1W	1207W25							
A&D ENGINEERING	1325A98	BALANCES	Analytical Balance 310gx0.1mg w/RS-232C	GR-300	7257-1W	1228T98									
A&D ENGINEERING	1325D45	BALANCES	Top Load Balance 210gx0.001g w/Breeze Break	GX-200	7258-1W	1228T99									
A&D ENGINEERING	1325D54	BALANCES	BALANCE 2100GX0.01G	GX-2000	7265-1W	1228U06									
A&D ENGINEERING	1325D48	BALANCES	Top Load Balance 410gx0.001g w/Breeze Break	GX-400	7258-1W	1228T99									
A&D ENGINEERING	1325D57	BALANCES	BALANCE 4100GX0.01G	GX-4000	7265-1W	1228U06									
A&D ENGINEERING	1325D63	BALANCES	BALANCE 6100GX0.1G	GX-6000		5kg	7010-1W	1207W21	200g	7016-1W	1207W25				
A&D ENGINEERING	1326A94	BALANCES	BALANCE 100G X .01G	HL-100		100g	7017-1W	0867B39	20g	7020-1W	0867B41				
A&D ENGINEERING	1226U09	BALANCES	SCALE IP65 3000g x 1g NTEP APPROVED	HL-3000WPN		2kg	7012-1W	1207W22							
A&D ENGINEERING	1325A15	BALANCES	BALANCE 120 X 0.0001G	HR-120	7255-1W	1228T96									
A&D ENGINEERING	1325A16	BALANCES	BALANCE 120X.0001G	HR-120-C	7255-1W	1228T96									
A&D ENGINEERING	1325A20	BALANCES	BALANCE 210 X 0.0001G	HR-200	7257-1W	1228T98									
A&D ENGINEERING	1325A23	BALANCES	BALANCE 210X.0001G	HR-200-C	7257-1W	1228T98									
A&D ENGINEERING	1325A10	BALANCES	Analytical Balance 60g x 0.1mg	HR-60	7252-1W	1228T93									
A&D ENGINEERING	1325A08	BALANCES	BALANCE 42/ 210X.00001/0.0001G	HR202C		200g	7016-1W	1207W25							
A&D ENGINEERING	1332B49	BALANCES	51X.002G MF50	MF-50		50g	7018-1W	0867B40	2g	7024-1W	0867B44				
A&D ENGINEERING	1332B46	BALANCES	51X.001G MX50	MX-50		50g	7018-1W	0867B40	2g	7024-1W	0867B44				
A&D ENGINEERING	1324Z56	BALANCES	POCKET BALANCE 200 X 0.1G	PV-200		200g	7016-1W	1207W25							
A&D ENGINEERING	1324Z59	BALANCES	POCKET BALANCE 500 X 0.1G	PV-500		500g	7014-1W	1207W24	200g	7016-1W	1207W25				
A&D ENGINEERING	1324Z50	BALANCES	POCKET BALANCE 60 X 0.01G	PV-60		50g	7018-1W	0867B40	20g	7020-1W	0867B41				
A&D ENGINEERING	1226U13	BALANCES	SCALE IP65 5000g x 2g NTEP APPROVED	SK-5000WP		5kg	7010-1W	1207W21							
A&D ENGINEERING	1326D01	BALANCES	BALANCE 120 X 0.01G	SV-120		100g	7017-1W	0867B39	20g	7020-1W	0867B41				
A&D ENGINEERING	1326D03	BALANCES	BALANCE 200 X 0.01G	SV-200		200g	7016-1W	1207W25	20g	7020-1W	0867B41				
A&D ENGINEERING	1326D07	BALANCES	BALANCE 610X0.1G	SV-610		500g	7014-1W	1207W24	200g	7016-1W	1207W25				
ADAM EQUIPMENT, INC	1224C91	BALANCES	CBK Bench Scale	CBK 16aH110v		5kg	7010-1W	1207W21	200g	7016-1W	1207W25				
ADAM EQUIPMENT, INC	1224C93	BALANCES	CBK Bench Scale	CBK 35a110v		10kg	7009-1W	1207W20	1kg	7013-1W	1207W23				
ADAM EQUIPMENT, INC	1224C90	BALANCES	CBK Bench Scale	CBK 8a110v		2kg	7012-1W	1207W22	200g	7016-1W	1207W25				
ADAM EQUIPMENT, INC	1224C98	BALANCES	CPWplus Bench Scale	CPWplus 35110v		20kg	7008-1W	1207W19							
ADAM EQUIPMENT, INC	1224D92	BALANCES	CPWplus P Bench Scale	CPWplus 35P220v		20kg	7008-1W	1207W19							
ADAM EQUIPMENT, INC	1224D93	BALANCES	CPWplus P Bench Scale	CPWplus 75P220v		50kg	7005-1W	1207W18							
ADAM EQUIPMENT, INC	1224C71	BALANCES	Core Compact Portable Balance	CQT 202110v		200g	7016-1W	1207W25	20g	7020-1W	0867B41				
ADAM EQUIPMENT, INC	1224C72	BALANCES	Core Compact Portable Balance	CQT 251110v		200g	7016-1W	1207W25							
ADAM EQUIPMENT, INC	1224D67	BALANCES	Dune Compact Balance	DCT 2000220v		2kg	7012-1W	1207W22							
ADAM EQUIPMENT, INC	1224C78	BALANCES	Dune Compact Balance	DCT 201110v		200g	7016-1W	1207W25							
ADAM EQUIPMENT, INC	1224C68	BALANCES	Highland Portable Precision Balance	HCB 1002110v	7263-1W	1228U04									

* If a Cal-Pak is not available we have provided recommended weights for balance linear verification.



Cal-Pak Reference Chart



800.345.2100 • 800.345.5232 fax
PO Box 99 • Swedesboro, NJ 08085
www.thomassci.com • value@thomassci.com



Vendor Name	Thomas part #	Category	Description	Manu Part #	Cal-Pak* #	Thomas	Recommended Weight #1			Recommended Weight #2			Recommended Weight #3		
							Size	Part #	Thomas	Size	Part #	Thomas	Size	Part #	Thomas
ADAM EQUIPMENT, INC	1224C63	BALANCES	Highland Portable Precision Balance	HCB 123110v	7256-1W	1228T97									
ADAM EQUIPMENT, INC	1224C65	BALANCES	Highland Portable Precision Balance	HCB 302110v			200g	7016-1W	1207W25						
ADAM EQUIPMENT, INC	1224C66	BALANCES	Highland Portable Precision Balance	HCB 602H110v			500g	7014-1W	1207W24	20g	7020-1W	0867B41			
ADAM EQUIPMENT, INC	1224C89	BALANCES	LBK Bench Scale	LBK 65a110v			20kg	7008-1W	1207W19						
ADAM EQUIPMENT, INC	1224C86	BALANCES	LBK Bench Scale	LBK 6a110v			2kg	7012-1W	1207W22						
ADAM EQUIPMENT, INC	1224D45	BALANCES	PGL Precision Balance	PGL 2002220v	7265-1W	1228U06									
ADAM EQUIPMENT, INC	1224C56	BALANCES	PGL Precision Balance	PGL 203110v	7258-1W	1228T99									
ADAM EQUIPMENT, INC	1219C33	BALANCES	PGL Precision Balance	PGL 3002	7265-1W	1228U06									
ADAM EQUIPMENT, INC	1219C34	BALANCES	PGL Precision Balance	PGL 303	7258-1W	1228T99									
ADAM EQUIPMENT, INC	1224C41	BALANCES	PGW Precision Balance	PGW 253e110v	7258-1W	1228T99									
ADAM EQUIPMENT, INC	1224C37	BALANCES	PW Analytical Balance	PW 124110v	7255-1W	1228T96									
ADAM EQUIPMENT, INC	1224C39	BALANCES	PW Analytical Balance	PW 254110v	7257-1W	1228T98									
ADAM EQUIPMENT, INC	1224C84	BALANCES	TBB Triple Beam Balance	TBB 2610S			500g	7014-1W	1207W24						
ADAM EQUIPMENT, INC	1224C83	BALANCES	TBB Triple Beam Balance	TBB 2610T			500g	7014-1W	1207W24						
ADAM EQUIPMENT, INC	1224D73	BALANCES	TBB Triple Beam Balance	TBB 610S			500g	7014-1W	1207W24						
ADAM EQUIPMENT, INC	1224C85	BALANCES	TBB Triple Beam Balance	TBB 610S			500g	7014-1W	1207W24						
ADAM EQUIPMENT, INC	1224D22	BALANCES	WBW Washdown Scale	WBW 9a110v			2kg	7012-1W	1207W22						
OHAUS SCALE CORPORATION	1368D00	BALANCES	BALANCE MODEL 1450-SD	1450-SD			2kg	7012-1W	1207W22	200g	7016-1W	1207W25			
OHAUS SCALE CORPORATION	1368D06	BALANCES	BALANCE MODEL 1550-SD	1550-SD			200g	7016-1W	1207W25						
OHAUS SCALE CORPORATION	1368D05	BALANCES	BALANCE MODEL 1560-SD	1560-SD			200g	7016-1W	1207W25						
OHAUS SCALE CORPORATION	1381H15	BALANCES	BALANCE MODEL 1650-00	1650-00			500g	7014-1W	1207W24						
OHAUS SCALE CORPORATION	1382F35	BALANCES	BALANCE MODEL 710-00	710-00			500g	7014-1W	1207W24						
OHAUS SCALE CORPORATION	1382F15	BALANCES	BALANCE MODEL 750-SO	750-SO			500g	7014-1W	1207W24						
OHAUS SCALE CORPORATION	1382F65	BALANCES	BALANCE MODEL 760-00	760-00			500g	7014-1W	1207W24						
OHAUS SCALE CORPORATION	1234X05	BALANCES	BALANCE ADVENT SL 3100GX0.01G	AS3102	7265-1W	1228U06									
OHAUS SCALE CORPORATION	1366T08	BALANCES	BALANCE ADVENTURER PRO 110GX0.1MG	AV114	7255-1W	1228T96									
OHAUS SCALE CORPORATION	1366T06	BALANCES	BALANCE ADVENTURER PRO 110GX0.1MG	AV114C	7255-1W	1228T96									
OHAUS SCALE CORPORATION	1366T36	BALANCES	BALANCE ADVENTURER PRO 2100GX0.1G	AV2101C			2kg	7012-1W	1207W22	200g	7016-1W	1207W25			
OHAUS SCALE CORPORATION	1366T28	BALANCES	BALANCE ADVENTURER PRO 2100GX0.01G	AV2102C	7265-1W	1228U06									
OHAUS SCALE CORPORATION	1366R80	BALANCES	BALANCE ADVENTURER PRO 210G X 0.001G	AV212			200g	7016-1W	1207W25	20g	7020-1W	0867B41			
OHAUS SCALE CORPORATION	1366T22	BALANCES	BALANCE ADVENTURER PRO 210GX0.01G	AV212C			200g	7016-1W	1207W25	20g	7020-1W	0867B41			
OHAUS SCALE CORPORATION	1366T16	BALANCES	BALANCE ADVENTURER PRO 210GX0.001G	AV213	7258-1W	1228T99									
OHAUS SCALE CORPORATION	1366T14	BALANCES	BALANCE ADVENTURER PRO 210GX0.001G	AV213C	7258-1W	1228T99									
OHAUS SCALE CORPORATION	1366T12	BALANCES	BALANCE ADVENTURER PRO 260GX0.1MG	AV264	7257-1W	1228T98									
OHAUS SCALE CORPORATION	1366T10	BALANCES	BALANCE ADVENTURER PRO 260GX0.1MG	AV264C	7257-1W	1228T98									
OHAUS SCALE CORPORATION	1366T32	BALANCES	BALANCE ADVENTURER PRO 3100GX0.01G	AV3102C	7265-1W	1228U06									
OHAUS SCALE CORPORATION	1366T18	BALANCES	BALANCE ADVENTURER PRO 310GX0.001G	AV313C	7258-1W	1228T99									
OHAUS SCALE CORPORATION	1366R88	BALANCES	BALANCE ADVENTURER PRO 4100G X 0.1G	AV4101			2kg	7012-1W	1207W22	200g	7016-1W	1207W25			
OHAUS SCALE CORPORATION	1366T38	BALANCES	BALANCE ADVEN PRO 4100GX0.1G	AV4101C			2kg	7012-1W	1207W22	200g	7016-1W	1207W25			
OHAUS SCALE CORPORATION	1366R82	BALANCES	BALANCE ADVENTURER PRO 410G X 0.01G	AV412			200g	7016-1W	1207W25	20g	7020-1W	0867B41			
OHAUS SCALE CORPORATION	1366T24	BALANCES	BALANCE ADVENTURER PRO 410GX0.01G	AV412C			200g	7016-1W	1207W25	20g	7020-1W	0867B41			
OHAUS SCALE CORPORATION	1366R78	BALANCES	BALANCE ADVENTURER PRO 53G X 0.001G	AV53			50g	7018-1W	0867B40	2g	7024-1W	0867B44			
OHAUS SCALE CORPORATION	1366T04	BALANCES	BALANCE ADVENTURER PRO 65GX0.1MG	AV64	7252-1W	1228T93									
OHAUS SCALE CORPORATION	1366R90	BALANCES	BALANCE ADVENTURER PRO 8100G X 0.1G	AV8101			5kg	7010-1W	1207W21	200g	7016-1W	1207W25			
OHAUS SCALE CORPORATION	1366T40	BALANCES	BALANCE ADVEN PRO 8100GX0.1G	AV8101C			5kg	7010-1W	1207W21	200g	7016-1W	1207W25			
OHAUS SCALE CORPORATION	1366R84	BALANCES	BALANCE ADVENTURER PRO 810G X 0.01G	AV812			500g	7014-1W	1207W24	20g	7020-1W	0867B41			
OHAUS SCALE CORPORATION	1366T26	BALANCES	BALANCE ADVENTURER PRO 810GX0.01G	AV812C			500g	7014-1W	1207W24	20g	7020-1W	0867B41			
OHAUS SCALE CORPORATION	1229X99	BALANCES	BALANCE CATAPULT 20 x 0.01kg	C11P20			20g	7020-1W	0867B41						
OHAUS SCALE CORPORATION	1229X96	BALANCES	BALANCE COMPACT 2000 x 1g	CL2000			2kg	7012-1W	1207W22						
OHAUS SCALE CORPORATION	1229X95	BALANCES	BALANCE COMPACT 200 x 0.1g	CL201			200g	7016-1W	1207W25						
OHAUS SCALE CORPORATION	1229X97	BALANCES	BALANCE COMPACT 5000 x 1g	CL5000			5kg	7010-1W	1207W21						
OHAUS SCALE CORPORATION	1367A71	BALANCES	BALANCE PRTBL 200G X 0.1G	CS200			200g	7016-1W	1207W25						
OHAUS SCALE CORPORATION	1367A76	BALANCES	BALANCE PRTBL 2000G X 1.0G	CS2000			2kg	7012-1W	1207W22						
OHAUS SCALE CORPORATION	1367A82	BALANCES	BALANCE PRTBL 5000G X 2.0G	CS5000			5kg	7010-1W	1207W21						
OHAUS SCALE CORPORATION	1229Y01	BALANCES	BALANCE DEF 3000 RECT BT 150 x 0.02kg	D31P150BX			100g	7017-1W	0867B39	20g	7020-1W	0867B41			
OHAUS SCALE CORPORATION	1229Y40	BALANCES	BALANCE D-5000 LP RT 100x0.02kg w/T51P	D51P100HL2			100g	7017-1W	0867B39	20g	7020-1W	0867B41			

* If a Cal-Pak is not available we have provided recommended weights for balance linear verification.



Cal-Pak Reference Chart



800.345.2100 • 800.345.5232 fax
PO Box 99 • Swedesboro, NJ 08085
www.thomassci.com • value@thomassci.com



Vendor Name	Thomas part #	Category	Description	Manu Part #	Cal-Pak* #	Thomas	Recommended Weight #1			Recommended Weight #2			Recommended Weight #3		
							Size	Part #	Thomas	Size	Part #	Thomas	Size	Part #	Thomas
OHAUS SCALE CORPORATION	1229Y42	BALANCES	BALANCE D-5000 LP RT 300x0.05kg w/T51P	D51P300HX2			200g	7016-1W	1207W25						
OHAUS SCALE CORPORATION	1232W82	BALANCES	BALANCE OHAUS DISCOVERY 110G	DV114C	7255-1W	1228T96									
OHAUS SCALE CORPORATION	1232W83	BALANCES	BALANCE OHAUS DISCOVERY 210G	DV214C	7257-1W	1228T98									
OHAUS SCALE CORPORATION	1232W84	BALANCES	BALANCE OHAUS DISCOVERY 81G	DV215CD	7257-1W	1228T98									
OHAUS SCALE CORPORATION	1366T07	BALANCES	BALANCE EXP. PRO 110G X 0.1MG W/O CAL	EP114	7255-1W	1228T96									
OHAUS SCALE CORPORATION	1366T05	BALANCES	BALANCE EXP. PRO W/INTERN-CALIB	EP114C	7255-1W	1228T96									
OHAUS SCALE CORPORATION	1366T65	BALANCES	BALANCE EXP. PRO 12000G X 0.1G W/CAL	EP12001C			10kg	7009-1W	1207W20	1kg	7013-1W	1207W23	200g	7016-1W	1207W25
OHAUS SCALE CORPORATION	1366T37	BALANCES	BALANCE 2100G X 0.01G W/CAL	EP2102C	7265-1W	1228U06									
OHAUS SCALE CORPORATION	1366T09	BALANCES	BALANCE ANALY W/CAL 210GX0.1MG	EP214C	7257-1W	1228T98									
OHAUS SCALE CORPORATION	1366T13	BALANCES	BALANCE EXP. PRO 100/210GX0.1/1MG W/CAL	EP214DC	7258-1W	1228T99									
OHAUS SCALE CORPORATION	1366T41	BALANCES	BALANCE 4100G X 0.01G	EP4102C	7265-1W	1228U06									
OHAUS SCALE CORPORATION	1366T21	BALANCES	BALANCE EXP. PRO 410G X 0.001G W/CAL	EP413C	7258-1W	1228T99									
OHAUS SCALE CORPORATION	1366T57	BALANCES	BALANCE EXP. PRO 6100G X 0.1G W/CAL	EP6101C			5kg	7010-1W	1207W21	200g	7016-1W	1207W25			
OHAUS SCALE CORPORATION	1366T47	BALANCES	BALANCE 6100G X 0.01G	EP6102	7267-1W	1228U08									
OHAUS SCALE CORPORATION	1366T45	BALANCES	BALANCE EXP. PRO 6100G X 0.01G W/CAL	EP6102C	7267-1W	1228U08									
OHAUS SCALE CORPORATION	1366T25	BALANCES	BALANCE EXP. PRO 610G X 0.001G W/CAL	EP613C	7261-1W	1228U02									
OHAUS SCALE CORPORATION	1367D62	BALANCES	BALANCE 66X0.05LB ES	ES30R			20kg	7008-1W	1207W19						
OHAUS SCALE CORPORATION	1218V32	BALANCES	Digital Balance, Precision, 2200g	EX2202	7265-1W	1228U06									
OHAUS SCALE CORPORATION	1218V26	BALANCES	Digital Balance, Analytical, 220g	EX224	7257-1W	1228T98									
OHAUS SCALE CORPORATION	1371J70	BALANCES	MOISTURE BALANCE MB45	MB45			50g	7018-1W	0867B40	2g	7024-1W	0867B44			
OHAUS SCALE CORPORATION	1229J32	BALANCES	BALANCE ANALYTICAL PIONEER 110x0.1mg	PA114	7255-1W	1228T96									
OHAUS SCALE CORPORATION	1229J37	BALANCES	BALANCE TOPLOADING PIONEER 1500x0.01g	PA1502	7263-1W	1228U04									
OHAUS SCALE CORPORATION	1229J34	BALANCES	BALANCE TOPLOADING PIONEER 150x0.001g	PA153	7256-1W	1228T97									
OHAUS SCALE CORPORATION	1229J33	BALANCES	BALANCE TOPLOADING PIONEER 210x0.1mg	PA214	7257-1W	1228T98									
OHAUS SCALE CORPORATION	1229J38	BALANCES	BALANCE TOPLOADING PIONEER 3100x0.01g	PA3102	7265-1W	1228U06									
OHAUS SCALE CORPORATION	1229J35	BALANCES	BALANCE TOPLOADING PIONEER 310x0.001g	PA313	7258-1W	1228T99									
OHAUS SCALE CORPORATION	1229J39	BALANCES	BALANCE TOPLOADING PIONEER 4100x0.1g	PA4101			2kg	7012-1W	1207W22	200g	7016-1W	1207W25			
OHAUS SCALE CORPORATION	1229J36	BALANCES	BALANCE TOPLOADING PIONEER 510x0.01g	PA512			50g	7018-1W	0867B40	20g	7020-1W	0867B41			
OHAUS SCALE CORPORATION	1229J31	BALANCES	BALANCE ANALYTICAL PIONEER 65x0.1mg	PA64	7252-1W	1228T93									
OHAUS SCALE CORPORATION	1367D20	BALANCES	BALANCE 3KGX0.2G	RD3RS			2kg	7012-1W	1207W22	200g	7016-1W	1207W25			
OHAUS SCALE CORPORATION	1367D74	BALANCES	BALANCE 6000X.1G RD6RM	RD6RM			5kg	7010-1W	1207W21						
OHAUS SCALE CORPORATION	1367D23	BALANCES	BALANCE 6KGX0.5G	RD6RS			5kg	7010-1W	1207W21	1kg	7013-1W	1207W23			
OHAUS SCALE CORPORATION	1380F12	BALANCES	SCHOOL BALANCE	SB1200			2kg	7012-1W	1207W22						
OHAUS SCALE CORPORATION	1367B18	BALANCES	BALANCE TOP LOAD 2000G X 0.1G	SP2001			2kg	7012-1W	1207W22	200g	7016-1W	1207W25			
OHAUS SCALE CORPORATION	1367B13	BALANCES	BALANCE TOP LOAD 200G X 0.01G	SP202			200g	7016-1W	1207W25	20g	7020-1W	0867B41			
OHAUS SCALE CORPORATION	1367B23	BALANCES	BALANCE TOP LOAD 4000G X 0.1G	SP4001			2kg	7012-1W	1207W22	200g	7016-1W	1207W25			
OHAUS SCALE CORPORATION	1367B16	BALANCES	BALANCE TOP LOAD 400G X 0.1G	SP401			200g	7016-1W	1207W25						
OHAUS SCALE CORPORATION	1367B14	BALANCES	BALANCE TOP LOAD 400G X 0.01G	SP402			200g	7016-1W	1207W25	20g	7020-1W	0867B41			
OHAUS SCALE CORPORATION	1367B19	BALANCES	BALANCE TOP LOAD 6000G X 1.0G	SP6000			5kg	7010-1W	1207W21						
OHAUS SCALE CORPORATION	1367B17	BALANCES	BALANCE TOP LOAD 600G X 0.1G	SP601			500g	7014-1W	1207W24						
OHAUS SCALE CORPORATION	1230J52	BALANCES	SCALE VALOR 1000 COMPACT 33LB	V11P15			10kg	7009-1W	1207W20						
OHAUS SCALE CORPORATION	1230J53	BALANCES	SCALE VALOR 1000 COMPACT 66LB	V11P30			20kg	7008-1W	1207W19						
OHAUS SCALE CORPORATION	1230J45	BALANCES	SCALE VALOR SS 3000G X 1G	V31X3			2kg	7012-1W	1207W22						
OHAUS SCALE CORPORATION	1230J43	BALANCES	SCALE VALOR SS 200G X .01G	V31XH202			200g	7016-1W	1207W25	20g	7020-1W	0867B41			
OHAUS SCALE CORPORATION	1366N61	BALANCES	BALANCE VOYAGER PRO W/CAL 110G	VP114CN	7255-1W	1228T96									
OHAUS SCALE CORPORATION	1366Q06	BALANCES	BALANCE VOYAGER PRO W/CAL 210G	VP213CN	7258-1W	1228T99									
OHAUS SCALE CORPORATION	1366N64	BALANCES	BALANCE VOYAGER PRO 100/210G	VP214DCN	7255-1W	1228T96									
OHAUS SCALE CORPORATION	1366Q14	BALANCES	BALANCE VOYAGER PRO CAL 4100	VP4102CN	7265-1W	1228U06									
SARTORIUS MECHANTRONICS	1229T73	BALANCES	ANALYTICAL BALANCE 110gx0.1mg	ALC-110.4	7255-1W	1228T96									
SARTORIUS MECHANTRONICS	1229T79	BALANCES	TOP LOADING BALANCE 2100x0.01g	ALC-2100.2	7265-1W	1228U06									
SARTORIUS MECHANTRONICS	1229T77	BALANCES	TOP LOADING BALANCE 810x0.01g	ALC-810.2			500g	7014-1W	1207W24	20g	7020-1W	0867B41			
SARTORIUS MECHANTRONICS	1229T93	BALANCES	TOPLOADING PORT SCALE 1500x0.1g	AY1501			1kg	7013-1W	1207W23	200g	7016-1W	1207W25			
SARTORIUS MECHANTRONICS	1229T88	BALANCES	TOPLOADING PORT SCALE 210x0.01g	AY212			200g	7016-1W	1207W25	20g	7020-1W	0867B41			
SARTORIUS MECHANTRONICS	1229T94	BALANCES	TOPLOADING PORT SCALE 3100x0.1g	AY3101			2kg	7012-1W	1207W22	200g	7016-1W	1207W25			
SARTORIUS MECHANTRONICS	1229T89	BALANCES	TOPLOADING PORT SCALE 410x0.01g	AY412			200g	7016-1W	1207W25	20g	7020-1W	0867B41			
SARTORIUS MECHANTRONICS	1229T95	BALANCES	TOPLOADING PORT SCALE 5100x0.1g	AY5101			5kg	7010-1W	1207W21	200g	7016-1W	1207W25			

* If a Cal-Pak is not available we have provided recommended weights for balance linear verification.



Cal-Pak Reference Chart



800.345.2100 • 800.345.5232 fax
PO Box 99 • Swedesboro, NJ 08085
www.thomassci.com • value@thomassci.com



Vendor Name	Thomas part #	Category	Description	Manu Part #	Cal-Pak* #	Thomas	Recommended Weight #1			Recommended Weight #2			Recommended Weight #3		
							Size	Part #	Thomas	Size	Part #	Thomas	Size	Part #	Thomas
SARTORIUS MECHANTRONICS	1229T91	BALANCES	TOPLOADING PORT SCALE 510x0.1g	AY511			500g	7014-1W	1207W24						
SARTORIUS MECHANTRONICS	1229T90	BALANCES	TOPLOADING PORT SCALE 610x0.01g	AY612			500g	7014-1W	1207W24	20g	7020-1W	0867B41			
SARTORIUS MECHANTRONICS	1352M64	BALANCES	BALANCE TE 120X0.1MG	AZ124	7255-1W	1228T96									
SARTORIUS MECHANTRONICS	1352M74	BALANCES	BALANCE TE 3100X0.01G	AZ1502	7265-1W	1228U06									
SARTORIUS MECHANTRONICS	1352M68	BALANCES	BALANCE TE 150X0.001G	AZ153	7256-1W	1228T97									
SARTORIUS MECHANTRONICS	1229T75	BALANCES	TOP LOADING BALANCE 150gx1mg	AZ153	7256-1W	1228T97									
SARTORIUS MECHANTRONICS	1352M79	BALANCES	BALANCE TE 210X0.01G	AZ212			200g	7016-1W	1207W25	20g	7020-1W	0867B41			
SARTORIUS MECHANTRONICS	1229T74	BALANCES	ANALYTICAL BALANCE 210gx0.1mg	AZ214	7257-1W	1228T98									
SARTORIUS MECHANTRONICS	1352M62	BALANCES	BALANCE TE 210X0.1MG	AZ214	7257-1W	1228T98									
SARTORIUS MECHANTRONICS	1352M71	BALANCES	BALANCE TE 3100X0.01G	AZ3102	7265-1W	1228U06									
SARTORIUS MECHANTRONICS	1352M65	BALANCES	BALANCE TE 310X0.001G	AZ313	7258-1W	1228T99									
SARTORIUS MECHANTRONICS	1229T76	BALANCES	TOP LOADING BALANCE 320gx1mg	AZ313	7258-1W	1228T99									
SARTORIUS MECHANTRONICS	1229T82	BALANCES	TOP LOADING BALANCE 4100x0.1g	AZ4101			2kg	7012-1W	1207W22	200g	7016-1W	1207W25			
SARTORIUS MECHANTRONICS	1352M85	BALANCES	BALANCE TE 610X0.1G	AZ601			500g	7014-1W	1207W24						
SARTORIUS MECHANTRONICS	1229T83	BALANCES	TOP LOADING BALANCE 6100x0.1g	AZ6101			5kg	7010-1W	1207W21	200g	7016-1W	1207W25			
SARTORIUS MECHANTRONICS	1352M91	BALANCES	BALANCE TE 640X.0001G	AZ64	7252-1W	1228T93									
SARTORIUS MECHANTRONICS	1352Q34	BALANCES	BALANCE .01G X 3200G	CP3202S	7265-1W	1228U06									
SARTORIUS MECHANTRONICS	1229U34	BALANCES	TOPLOADER BALANCE 10kg x 0.1g	CPA10001			10kg	7009-1W	1207W20	1kg	7013-1W	1207W23	200g	7016-1W	1207W25
SARTORIUS MECHANTRONICS	1229U16	BALANCES	ANALYTICAL BALANCE 120g x 0.1mg	CPA124S	7255-1W	1228T96									
SARTORIUS MECHANTRONICS	1229U25	BALANCES	TOPLOADER BALANCE 2200g x 0.01g	CPA2202S	7265-1W	1228U06									
SARTORIUS MECHANTRONICS	1229U19	BALANCES	MILLIGRAM BALANCE 220g x 0.001g	CPA223S	7258-1W	1228T99									
SARTORIUS MECHANTRONICS	1229U17	BALANCES	ANALYTICAL BALANCE 220g x 0.1mg	CPA224S	7257-1W	1228T98									
SARTORIUS MECHANTRONICS	1229U14	BALANCES	SEMIMICROBALANCE 100/200g x .01/1mg	CPA225D	7254-1W	1228T95									
SARTORIUS MECHANTRONICS	1229U12	BALANCES	MICROBALANCE 0.51/2g x .001/.002/.005mg	CPA2P			500mg	7026-1W	0867B46	50mg	7030-1W	0867B49	2mg	7036-1W	0867B53
SARTORIUS MECHANTRONICS	1229U26	BALANCES	TOPLOADER BALANCE 3200g x 0.01g	CPA3202S	7265-1W	1228U06									
SARTORIUS MECHANTRONICS	1229U20	BALANCES	MILLIGRAM BALANCE 320g x 0.001g	CPA323S	7258-1W	1228T99									
SARTORIUS MECHANTRONICS	1229U18	BALANCES	ANALYTICAL BALANCE 320g x 0.1mg	CPA324S	7257-1W	1228T98									
SARTORIUS MECHANTRONICS	1229U37	BALANCES	SCALE 8/16/34kg x 0.1/0.2/0.5g	CPA34001P			5kg	7010-1W	1207W21	200g	7016-1W	1207W25			
SARTORIUS MECHANTRONICS	1229U38	BALANCES	HIGH CAPACITY SCALE 34kg x 0.1g	CPA34001S			20kg	7008-1W	1207W19	2kg	7012-1W	1207W22	200g	7016-1W	1207W25
SARTORIUS MECHANTRONICS	1229U27	BALANCES	TOPLOADER BALANCE 4200g x 0.01g	CPA4202S	7265-1W	1228U06									
SARTORIUS MECHANTRONICS	1229U15	BALANCES	ANALYTICAL BALANCE 64g x 0.1mg	CPA64	7252-1W	1228T93									
SARTORIUS MECHANTRONICS	1229U01	BALANCES	PORTABLE SCALE 2100x1g	EC-2100			2kg	7012-1W	1207W22						
SARTORIUS MECHANTRONICS	1229T99	BALANCES	PORTABLE SCALE 210x0.1g	EC-211			200g	7016-1W	1207W25						
SARTORIUS MECHANTRONICS	1229U02	BALANCES	PORTABLE SCALE 4100x1g	EC-4100			2kg	7012-1W	1207W22						
SARTORIUS MECHANTRONICS	1229U00	BALANCES	PORTABLE SCALE 410x0.1g	EC-411			200g	7016-1W	1207W25						
SARTORIUS MECHANTRONICS	1229T64	BALANCES	PREC BAL w/MOTOR INT CALIB 120gx0.1mg	ED124S	7255-1W	1228T96									
SARTORIUS MECHANTRONICS	1229T65	BALANCES	PREC BAL w/MOTOR INT CALIB 220gx0.1mg	ED224S	7257-1W	1228T98									
SARTORIUS MECHANTRONICS	1353K67	BALANCES	BALANCE 2200 X 0.01 G	FC2CCE-SX	7265-1W	1228U06									
SARTORIUS MECHANTRONICS	1229U07	BALANCES	INDUSTRIAL PRECISION SCALE 20kpx2g	IW2P1U-30ED-L			20kg	7008-1W	1207W19						
SARTORIUS MECHANTRONICS	1362A60	BALANCES	MICROBALANCE 5.1G X 1.0UG	ME5			5g	7022-1W	0867B43	500mg	7026-1W	0867B46	2mg	7036-1W	0867B53
SARTORIUS MECHANTRONICS	1229U03	BALANCES	POCKET PRO SCALE 210x0.1g	PP-201			200g	7016-1W	1207W25						
SARTORIUS MECHANTRONICS	1229U04	BALANCES	POCKET PRO SCALE 410x0.01g	PP-401			200g	7016-1W	1207W25	20g	7020-1W	0867B41			
SARTORIUS MECHANTRONICS	0811B37	BALANCES	BALANCE ULTRAMICRO 2.1G X0.1UG	SE2			2g	7024-1W	0867B44						
SARTORIUS MECHANTRONICS	1229T85	BALANCES	MILLIGRAM SCALE 300gx1mg	VIC-303	7258-1W	1228T99									
SARTORIUS STEDIM	1229T72	BALANCES	ANALYTICAL BALANCE 80gx0.1mg	ALC-80.4	7252-1W	1228T93									
SARTORIUS STEDIM	1352M76	BALANCES	BALANCE TE 610X0.01G	TE612			500g	7014-1W	1207W24	20g	7020-1W	0867B41			
Thomas Brand	1233Z73	BALANCES	BALANCE 102G X 0.001G	165-9327-310	7256-1W	1228T97									
Thomas Brand	1233Z74	BALANCES	BALANCE 610G X 0.01G	165-9423-310			500g	7014-1W	1207W24	20g	7020-1W	0867B41			
Thomas Brand	1233Z75	BALANCES	BALANCE 1020G X 0.01G	165-9425-310	7263-1W	1228U04									
Thomas Brand	1233Z76	BALANCES	BALANCE 4100G X 0.1G	165-9514-310			2kg	7012-1W	1207W22	200g	7016-1W	1207W25			
Thomas Brand	1233Z77	BALANCES	BALANCE 8100G X 0.1G	165-9518-310			5kg	7010-1W	1207W21	200g	7016-1W	1207W25			
Thomas Brand	1233Z78	BALANCES	BALANCE 12100G X 1.0G	165-9604-310			10kg	7009-1W	1207W20	1kg	7013-1W	1207W23	200mg	7028-1W	0867B47
Thomas Brand	1324Y02	BALANCES	BALANCE 120G X 0.1MG	320-9200-317	7255-1W	1228T96									
Thomas Brand	1324Y04	BALANCES	BALANCE 220G X 0.1MG	320-9201-317	7257-1W	1228T98									
Thomas Brand	1324Y08	BALANCES	BALANCE 120G X 0.1MG	320-9203-313	7255-1W	1228T96									
Thomas Brand	1324Y10	BALANCES	BALANCE 220G X 0.1MG	320-9204-313	7257-1W	1228T98									
Thomas Brand	1324Y22	BALANCES	BALANCE 320G X 1MG	320-9307-313	7258-1W	1228T99									

* If a Cal-Pak is not available we have provided recommended weights for balance linear verification.



Cal-Pak Reference Chart



800.345.2100 • 800.345.5232 fax
PO Box 99 • Swedesboro, NJ 08085
www.thomassci.com • value@thomassci.com



Vendor Name	Thomas part #	Category	Description	Manu Part #	Cal-Pak* #	Thomas	Recommended Weight #1			Recommended Weight #2			Recommended Weight #3		
							Size	Part #	Thomas	Size	Part #	Thomas	Size	Part #	Thomas
Thomas Brand	1324Y16	BALANCES	BALANCE 2200G X 10MG	320-9401-317	7265-1W	1228U06									
Thomas Brand	1324Y26	BALANCES	BALANCE 2200G X 10MG	320-9409-313	7265-1W	1228U06									
Thomas Brand	1324Y28	BALANCES	BALANCE 4200G X 10MG	320-9410-313	7265-1W	1228U06									
Thomas Brand	1324Y20	BALANCES	BALANCE 6200G X 0.1G	320-9501-317			5kg	7010-1W	1207W21	200g	7016-1W	1207W25			
Thomas Brand	1324Y30	BALANCES	BALANCE 6200G X 10MG	320-9507-313	7267-1W	1228U08									
Thomas Brand	0782A02	BALANCES	TXS BAL 220X.1MG 230V	3209204323	7257-1W	1228T98									

* If a Cal-Pak is not available we have provided recommended weights for balance linear verification.