


GEB Universal tips for TWC thin wallcrown And for improved flexibility.

Product details at-a-glance

GEB Cat#	Max Volume	Tip Color	Low Retention	Sterilized	Filtered (2.5*3.5)mm	Packaging Spec.
FT0020-B-N	20ul	Natural			Filtered	1000tips/bag,
FT0020-B-N-LB	20ul	Natural	Low Retention		Filtered	10bags/case, 10000tips/case
Case dimension L*W*H: 450*265*270mm.			Gross Weight: 4.7kg.			
FT0020-R-N	20ul	Natural			Filtered	96tips/rack,
FT0020-R-N-LB	20ul	Natural	Low Retention		Filtered	10racks/pack,
FT0020-R-NS	20ul	Natural		Sterilized	Filtered	5packs/case,
FT0020-R-NS-LB	20ul	Natural	Low Retention	Sterilized	Filtered	4800tips/case.
Case dimension L*W*H: 500*265*450mm.			Gross Weight: 7.8 kg.			

- ✧ 20μL with graduation marks at 10ul
- ✧ 96 Tip/Hanging Tip Rack format / Bulk
- ✧ Polypropylene
- ✧ Presterile/Nonsterile
- ✧ Filtered (With aerosol-resistant filters) /NonFiltered
- ✧ Certified RNase-/DNase-free
- ✧ Nonpyrogenic
- ✧ Endotoxin-Safe meets USP/FDA standards
- ✧ Comply With USP 3127 (Heavy Metal Contamination) standards
- ✧ Compatibility: Universal tips workstations.

Drawing:



Prepared by: Anne Oct.19, 2020 Approve: WeiMing Yu Oct.19,2020