




Eisco Industrial offers a vast selection of premium industrial quality labware. Our products have been carefully designed by laboratory professionals along-side experienced engineers with quality and durability in mind. They are manufactured with only the highest-grade raw materials and according to international quality standards. From large capacity glassware to sturdy clamps and lab stands you will find the high quality products that you need for botanical extraction laboratories of all sizes.

### Conical Filtering Flasks

Heavy wall conical filtering flasks with integral glass side arm. Made from high quality borosilicate 3.3 glass. Autoclavable.



### Conical Filtering Flasks


Thomas No.	Mfr. No.	Capacity	Graduations
23A00Y642	CH0420C	500 ml	100 ml
23A00Y643	CH0420D	1000 ml	200 ml
23A00Y644	CH0420E	2000 ml	250 ml
23A00Y645	CH0420F	5000 ml	500 ml

Additional product sizes and styles are available.

Please visit [ThomasSci.com](http://ThomasSci.com) for additional products and product specifications.

## Erlenmeyer Flasks


Narrow neck Erlenmeyer flasks made from high quality borosilicate 3.3 glass. DIN ISO 1773 compliant. Autoclavable.



Erlenmeyer Flasks			
Thomas No.	Mfr. No.	Capacity	Graduations
24A00A551	CH0424G	500 ml	100 ml
24A00A552	CH0424H	1000 ml	200 ml
24A00A553	CH0424I	2000 ml	250 ml
24A00A554	CH0424J	3000 ml	500 ml
24A00A555	CH0424K	5000 ml	500 ml

## Low Form Beakers

Low form beakers made from high quality borosilicate 3.3 glass. Graduated as per DIN 12231 ISO 3819. Autoclavable.



Low Form Beakers			
Thomas No.	Mfr. No.	Capacity	Graduations
23A00X862	CH0126I	500 ml	50 ml
23A00X863	CH0126J	600 ml	50 ml
23A00X864	CH0126J8	800 ml	100 ml
23A00X865	CH0126K	1000 ml	100 ml
23A00X866	CH0126L	2000 ml	100 ml
23A00X867	CH0126M	3000 ml	250 ml
23A00X868	CH0126N	5000 ml	500 ml
23A00X869	CH0126O	10,000 ml	1,000 ml

## Reagent Bottles

Borosilicate 3.3 glass reagent bottles with screw caps (GL45) and pouring rings. Both bottle and cap are autoclavable. DIN ISO 4796-1 compliant.



Reagent Bottles				
Clear		Amber		Capacity
Thomas No.	Mfr. No.	Thomas No.	Mfr. No.	
24A00A559	CH0164A	24A00B241	CH0165A	100 ml
24A00A560	CH0164B	24A00B242	CH0165B	250 ml
24A00A561	CH0164C	24A00B243	CH0165C	500 ml
24A00A562	CH0164D	24A00B244	CH0165D	1000 ml
24A00A563	CH0164E	24A00B245	CH0165E	2000 ml

Additional product sizes and styles are available.

Please visit [ThomasSci.com](https://www.thomasci.com) for additional products and product specifications.

## Squibb Separatory Funnels

Squibb separatory funnels with interchangeable polypropylene stoppers and key stopcocks. Made from high quality borosilicate 3.3 glass.

Squibb Separatory Funnels					
Glass Stopcock		PTFE Stopcock		Capacity	Socket Size
Thomas No.	Mfr. No.	Thomas No.	Mfr. No.		
24A00A564	CH0476C	23A00Y065	CH0479D	500 ml	24/29
24A00A565	CH0476D	23A00Y066	CH0479E	1000 ml	29/32
24A00A566	CH0476E	23A00Y068	CH0479F	2000 ml	29/32



## Flat Bottom Boiling Flasks

Flat bottom boiling flasks made from high quality borosilicate 3.3 glass. Narrow neck design with a beaded rim for increased durability. DIN ISO 1773 compliant. Autoclavable.



Flat Bottom Boiling Flasks		
Thomas No.	Mfr. No.	Capacity
23A00Y633	CH0410F	500 ml
23A00Y634	CH0410G	1000 ml
23A00Y635	CH0410H	2000 ml
24A00A567	CH0410I	3000 ml
23A00Y636	CH0410J	5000 ml
23A00Y637	CH0410K	6000 ml
23A00Y638	CH0410L	10,000 ml
23A00Y639	CH0410M	20,000 ml

## 3 Neck Distilling Flasks

3 Neck distilling flasks made from high quality borosilicate 3.3 glass feature round bottoms, ground joints and 2 oblique necks. Autoclavable.



3 Neck Distilling Flasks				
Thomas No.	Mfr. No.	Capacity	Center Socket	Side Socket
24A00A568	CH0418P	500 ml	24/29	14/23
24A00A569	CH0418Q	500 ml	24/29	19/26
24A00A570	CH0418R	500 ml	29/32	19/26
24A00A571	CH0418S	1000 ml	24/29	14/23
24A00A572	CH0418T	1000 ml	24/29	19/26
24A00A573	CH0418U	1000 ml	29/32	19/26
24A00A574	CH0418V	2000 ml	24/29	19/26
24A00A575	CH0418W	2000 ml	29/32	24/29

Additional product sizes and styles are available.

Please visit [ThomasSci.com](https://www.ThomasSci.com) for additional products and product specifications.

## Stainless Steel Laboratory Jacks

Premium laboratory jacks feature stainless steel top plates, slip proof rubber feet, and a screw operated mechanism that is easy to adjust.

Stainless Steel Laboratory Jacks					
Thomas No.	Mfr. No.	Top Plate	Max. Open Ht.	Dynamic Load	Static Load
23A00W301	SSJCT30	8 x 8 in	11 in	15.4 lbs	66 lbs
23A00W302	SSJCT40	10 x 10 in	13 in	15.4 lbs	66 lbs
23A00W303	SSJCT50	12 x 12 in	19 in	15.4 lbs	66 lbs
23A00W304	SSJCT60	16 x 16 in	20 in	15.4 lbs	66 lbs



## Premium Laboratory Clamps

Premium lab clamps feature solid zinc die-cast construction with Zamac 3 Alloy + aluminum alloy. They are coated with a high grade, non-toxic powder finish, are rust resistant and durable. Please visit [thomassci.com](http://thomassci.com) for a complete selection of clamps and specifications.



Premium Laboratory Clamps		
Thomas No.	Mfr. No.	Capacity
23A00W318	PRCMP01	Premium Clamp, 3 Finger, Vinyl Coated
23A00W319	PRCMP01BH	Premium Clamp, 3 Finger, Vinyl Coated, W/ Bosshead
23A00W320	PRCMP02	Premium Clamp, 4 Prong, Cork Lined
23A00W321	PRCMP02BH	Premium Clamp, 4 Prong, Cork Lined, W/ Bosshead
23A00W322	PRCMP03	Premium Clamp, 4 Prong, Vinyl Coated
23A00W323	PRCMP03BH	Premium Clamp, 4 Prong, Vinyl Coated, W/ Bosshead
23A00W324	PRCMP04	Premium Clamp, 2 Prong, Vinyl Coated, 2x Adj
23A00W325	PRCMP04BH	Premium Clamp, 2 Prong, Vinyl Coated, W/ Bosshead, 2x Adj
23A00W326	PRCMP05	Premium Clamp, 2 Prong, Vinyl Coated
23A00W327	PRCMP05BH	Premium Clamp, 2 Prong, Vinyl Coated, W/ Bosshead
23A00W328	PRCMP06	Premium Clamp, 3 Prong, Vinyl Coated, 2x Adj
23A00W329	PRCMP06BH	Premium Clamp, 3 Prong, Vinyl Coated, W/Bosshead, 2x Adj
23A00W330	PRCMP07	Premium Clamp, Long, 3 Finger, Vinyl Coated
23A00W331	PRCMP07BH	Premium Clamp, Long, 3 Finger, Vinyl Coated, W/ Bosshead
23A00W332	PRCMP08	Premium Clamp, 2 Prong, Cork Lined
23A00W333	PRCMP08BH	Premium Clamp, 2 Prong, Cork Lined, W/ Bosshead
23A00W334	PRCMP09	Premium Clamp, 3 Finger, Cork Lined
23A00W335	PRCMP09BH	Premium Clamp, 3 Finger, Cork Lined, W/ Bosshead

## Heavy Duty Retort Base and Rod Sets

These sets have been created with premium quality metalware and will provide maximum stability and durability. The bases are cast iron with a high grade non-toxic powder coating. The stainless steel rods have precise threading that results in a snug and stable fit into the bases.

Heavy Duty Retort Base and Rod Sets			
Thomas No.	Mfr. No.	Base	Rod
23A00W310	CH0653D1RD3	6 x 9 in	24 in (0.4 in dia)
23A00W311	CH0653E1RD4	6 x 11 in	35.5 in (0.5 in dia)



Additional product sizes and styles are available.

Please visit [ThomasSci.com](http://ThomasSci.com) for additional products and product specifications.