

PMM1 Polyclonal Antibody

Catalog number: 11252-2-AP


Size: 25 µg/150 µl

Source: Rabbit

Isotype: IgG

Synonyms:

PMM1; Sec53



human brain tissue were subjected to SDS PAGE followed by western blot with 11252-2-AP(PMM1 Antibody) at dilution of 1:400

Applications

Tested applications:	ELISA, WB
Cited applications:	WB
Species specificity:	Human, mouse, rat; other species not tested.
Cited species:	Mouse
Positive WB detected in:	Human brain tissue
Calculated PMM1 MW:	29kd
Observed PMM1 MW:	29kd
Recommended dilution:	WB: 1:200-1:1000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation, FC = Flow cytometry

Immunogen information

Immunogen:	Ag1782
GenBank accession number:	BC016818
Gene ID (NCBI):	5372
Full name:	Phosphomannomutase 1

Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.1% sodium azide and 50% glycerol pH 7.3. Store at -20°C.