


Description: Face Mask with Earloops (Model No. XC3000-5B)
Size: 17.5*9.5cm
Packaging: 50/PK; 10 PKS/CS; 500/CS
Color: Blue
Material: 25g PP + 25g filter paper + 25g PP
Specifications: Bacterial Filtration Efficiency (BFE): 98.5%
Designed for use in ISO 6 cleanroom environments



Item	Size	Tolerance	Comment
A	17.5cm	+/-1cm	
B	9.5cm	+/-1cm	
C	11cm	+/-1cm	
D	16cm	+/-2cm	

Note: All the unmarked tolerance is ± 2 cm.

Test Item	Result	Test Standard
Bacterial Filtration Efficiency (BFE)	>99%	ASTM F2101
Differential Pressure (mm H ₂ O/cm ²)	5.28	
Latex Particle Challenge (PFE)	>98.6%	ASTM F2299/F2299M-03