

## NMNAT3 Polyclonal Antibody

Catalog number: 13236-1-AP

Size: 37 µg/150 µl

Source: Rabbit


Isotype: IgG

Synonyms:

NMNAT3; NaMN

adenylyltransferase 3,

NMNAT3, PNAT 3, PNAT3



mouse kidney tissue were subjected to SDS PAGE followed by western blot with 13236-1-AP (NMNAT3 Antibody) at dilution of 1:300

### Applications

Tested applications:	ELISA, WB
Species specificity:	Human, Mouse, Rat; other species not tested.
Calculated NMNAT3 MW:	215aa, 24 kDa
Observed NMNAT3 MW:	24-28 kDa
Positive WB detected in	Mouse kidney tissue, mouse skeletal muscle tissue
Recommended dilution:	WB: 1:200-1:1000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation

### Immunogen information

Immunogen:	Ag3965
GenBank accession number:	BC034374
Gene ID (NCBI):	349565
Full name:	Nicotinamide nucleotide adenylyltransferase 3

### Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.1% sodium azide and 50% glycerol pH 7.3. Store at -20°C.