

SMA4 Polyclonal Antibody

Catalog number: 10318-1-AP

Size: 20 µg/150 µl

Source: Rabbit


Isotype: IgG

Synonyms:

SMA4; FLJ36702; FLJ39772;

MGC22265; MGC60382; SMA3;

b55C20.2



HeLa cells were subjected to SDS PAGE followed by western blot with 10318-1-AP(SMA4 Antibody) at dilution of 1:300

Applications

Tested applications:	ELISA, WB
Species specificity:	Human; other species not tested.
Positive WB detected in:	HeLa cells
Calculated SMA4 MW:	15kd
Observed SMA4 MW:	15kd
Recommended dilution:	WB: 1:200-1:1000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation, FC = Flow cytometry

Immunogen information

Immunogen:	Ag0372
GenBank accession number:	BC002622
Gene ID (NCBI):	11039
Full name:	Glucuronidase, beta pseudogene

Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.1% sodium azide and 50% glycerol pH 7.3. Store at -20°C.