



Cryovials

For over 25 years Crystal Technology & Industries has been driven to deliver exceptional value to the scientific marketplace. Our engineers collaborate with top-tier research institutions to design & develop innovative Benchtop Equipment and Cryo Storage Solutions for today's Scientist.

Crystal Clear Leak Proof Cryovials for Biological Sample Storage at Temperatures Down to -196°C

IN STOCK NOW



Product Features:

- Caps and Vials are molded from Medical Grade Polypropylene ensuring the optimum seal to prevent leakage
- The Internal Thread Screw Cap design & O-ring reduces contamination
- Sterile (Gamma Radiation) and free of endotoxin, RNase, DNase, DNA, PCR inhibitors, and pyrogens
- Rounded internal bottom of vial promotes complete pour out
- Crystal Clear Transparent vial for easy observation
- Large white writing area allows easy marking
- Manufactured in Cleanroom environment
- Designed for freezing samples in (vapor phase) liquid nitrogen
Not to be used in the liquid phase of liquid nitrogen

Specifications:

Part Number	Material	Outer Dimension	Volume Capacity	Temperature Range	Qty. in PK	Description
CVI-2.0-B	PP	Ø 12.5 x 50mm	2.0ml	-196 ~ 121°C	50	Internal Thread, Self Standing