

Product Name:	ABS LD 20
Description:	ABS LD Class 20 Liter Liquid Nitrogen Capacity 0.18 Liter/Day Static Evaporation Rate Neck Diameter (mm): 56 Exterior Height (mm): 622.3 Exterior Diameter: 368.3 Designed for storage of LN2 Low LN2 consumption Durable composite lid High strength neck tube reduces LN2 loss Strong, lightweight aluminum construction Safer operation & maintenance One-year parts & labor warranty Three-year vacuum warranty

