

TNFSF12 Polyclonal Antibody

Catalog number: 12537-1-AP

Size: 44 µg/150 µl

Source: Rabbit

Isotype: IgG

Synonyms:

TNFSF12; APO3 ligand, APO3L,
DR3LG, TNFSF12, TWEAK

Applications

Tested applications: ELISA
Species specificity: Human; other species not tested.
Calculated TNFSF12 MW: 330aa,37 kDa
Observed TNFSF12 MW:

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation

Immunogen information

Immunogen: Ag3242
GenBank accession number: BC019047
Gene ID (NCBI): 8742
Full name: Tumor necrosis factor (ligand) superfamily, member 12

Product information

Purification method: Antigen affinity purification
Storage: PBS with 0.1% sodium azide and 50% glycerol pH 7.3. Store at -20°C.