

## GLYAT Polyclonal Antibody

Catalog number: 10900-1-AP

Size: 22 µg/150 µl

Source: Rabbit


Isotype: IgG

Synonyms:

GLYAT; AAc, ACGNAT, CAT, GAT,

GLYAT, glycine N

acyltransferase, HRP 1(CLP)



mouse liver tissue were subjected to SDS PAGE followed by western blot with 10900-1-AP( GLYAT Antibody) at dilution of 1:600

### Applications

Tested applications:	ELISA, WB
Species specificity:	Human, Mouse; other species not tested.
Calculated GLYAT MW:	30 kDa
Observed GLYAT MW:	35 kDa
Positive WB detected in	Mouse liver tissue
Recommended dilution:	WB: 1:200-1:2000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation

### Immunogen information

Immunogen:	Ag1351
GenBank accession number:	BC009785
Gene ID (NCBI):	10249
Full name:	Glycine-N-acyltransferase

### Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.1% sodium azide and 50% glycerol pH 7.3. Store at -20°C.