

## B4GALNT1 Polyclonal Antibody

Catalog number: 13396-1-AP


Size: 49 µg/150 µl

Source: Rabbit

Isotype: IgG

Synonyms:

B4GALNT1; B4GALNT1,  
GALGT, GalNAc T, GALNACT,  
GM2/GD2 synthase, SIAT2



COLO 320 cells were subjected to SDS PAGE followed by western blot with 13396-1-AP(B4GALNT1 antibody) at dilution of 1:400

### Applications

Tested applications:	ELISA, WB
Species specificity:	Human; other species not tested.
Calculated B4GALNT1 MW:	532aa,58 kDa
Observed B4GALNT1 MW:	67.5kd
Positive WB detected in	COLO 320 cells, A375 cells
Recommended dilution:	WB: 1:500-1:5000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation

### Immunogen information

Immunogen:	Ag4065
GenBank accession number:	BC029828
Gene ID (NCBI):	2583
Full name:	Beta-1,4-N-acetyl-galactosaminyl transferase 1

### Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3. Store at -20°C.