

## Calpain 7 Polyclonal Antibody

Catalog number: 13870-1-AP

Size: 20 µg/150 µl

Source: Rabbit


Isotype: IgG

Synonyms:

CAPN7; calpain 7, CALPAIN7,

CAPN7, FLJ36423, PalB

homolog, PALBH



human heart tissue were subjected to SDS PAGE followed by western blot with 13870-1-AP(CAPN7 antibody) at dilution of 1:1000

### Applications

|                          |  |
|--------------------------|--|
| Tested applications:     | ELISA, WB                                  |
| Species specificity:     | Human,Mouse,Rat; other species not tested. |
| Calculated Calpain 7 MW: | 93kd                                       |
| Observed Calpain 7 MW:   | 92 kDa                                     |
| Positive WB detected in  | Human heart tissue                         |
| Recommended dilution:    | WB: 1:200-1:2000                           |

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation

### Immunogen information

|                           |           |
|---------------------------|-----------|
| Immunogen:                | Ag4820    |
| GenBank accession number: | BC056202  |
| Gene ID (NCBI):           | 23473     |
| Full name:                | Calpain 7 |

### Product information

|                      |  |
|----------------------|--|
| Purification method: | Antigen affinity purification  |
| Storage:             | PBS with 0.02% sodium azide and 50% glycerol pH 7.3. Store at -20°C. |