


The Multi-Task Genie™ Family

THE ORIGINAL

TIMER

DIGITAL



U.S. Pats. 4,781,487 / 5,707,861

SPECIFICATIONS

Vortex-Genie Model	Vortex-Genie® 2	Vortex-Genie® 2T	Digital Vortex-Genie® 2
Supplied	3-inch Platform, Pop-off Cup	3-inch Platform, Pop-off Cup	3-inch Platform, Pop-off Cup
Speed (RPM)	600-3200	600-3200	500-3000
Timer - Touch On: Hands-Free:	N/A	1-60 seconds 1-60 minutes or continuous	1-99 seconds 1-99 minutes or continuous
Weight	8.8 lb	8.8 lb	8.8 lb
Dimensions (base) (D x W x H)	6.5 x 4.8 x 6.5in	6.5 x 4.8 x 6.5in	6.5 x 4.8 x 6.5in
Markings / Warranty	ETL, CE / 2 years	ETL, CE / 2 years	CE / 2 years
CATALOG NUMBERS			
120V (North America)	8294D17	8294D22	1226W53

Multi-Task Genie™ Family Accessories



Fig.	Cat. No.	Description	Fig.	Cat. No.	Description
2,3	8294D35	Model H301 Multiple Sample Starter Set	13	1220A01	TurboMix™ Attachment
--	8294D65	6-inch Platform	14	1220A00	Universal Microtube Holder - (Requires 1220A01)
2	8294D80	Microtube Foam Insert (2) - (Requires 8294D65)	15	8294D82	Large Ampule/Tube Attachment
3	8294D70	Microplate Foam Insert (2) - (Requires 8294D65)	16	8294D81	Small Ampule/Tube Attachment
5-10	8294D37	Model H302 Large Sample Set	17	1224N32	Horizontal 50ml Tube Holder (6 tubes)
5	8294D48	Recessed Platform without Bands	18	1224N33	Horizontal 15ml Tube Holder (12 tubes)
6	8294D71	9-16mm Tube Foam Inserts (2) - (Requires 8294D48)	19	1224N34	Horizontal Microtube Holder (24 tubes)
7	8294D72	14-29mm Tube Foam Inserts (2) - (Requires 8294D48)	20	1216T73	Vertical High Speed 50ml Tube Holder (3 tubes)
8	8294D73	29-37mm Tube Foam Inserts (2) - (Requires 8294D48)	21	1216T74	Vertical High Speed 15ml Tube Holder (6 tubes)
9	8294D74	Blank Foam Inserts (2) - (Requires 8294D48)	22	1216T75	Vertical High Speed Microtube Holder (9 tubes)
6-9	8294D76	Assorted Foam Inserts (4) - (Requires 8294D48)	23	1224N35	Vertical 50ml Tube Holder (6 tubes)
10	8294D69	Elastic Bands (2) - (Requires 8294D48)	24	1224N36	Vertical 15ml Tube Holder (12 tubes)
11	8294D40	Pop-off Cup	25	1224N37	Vertical Microtube Holder (24 tubes)
12	8294D77	3-inch Platform with Rubber Cover	26	1221S04	Vertical (6) 50ml, (6) 15ml, (12) Microtube Holder
--	8294D47	Rubber Cover only for 3-inch Platform	--	1218A76	Bottom Closure with Screws Kit - (Inquire for older models)



800.345.2100 • 800.345.5232 fax
 PO Box 99 • Swedesboro, NJ 08085
 www.thomassci.com • value@thomassci.com

