

Best Seller!

KingFisher™ Compatible Deep Well Plates


- Printed black alpha-numerics
- Square wells and conical v-bottoms
- Available sterile or non-sterile

labForce® multiwell plates are fully compatible with KingFisher™ automation platforms, as well as Accuris IsoPure™. All four plates feature square wells with conical v-bottoms and printed alpha-numerics. They are available in sterile or non-sterile configurations.



23A00W712
96×2.2mL wells

Unique Feature!
Printed Alpha-Numerics



23A00W668
24×10mL wells



23A00W682
96×0.5mL wells



23A00W692
96×1.0mL wells

KingFisher™ Compatible Deep Well Plates - Sterile

Cat. No.	Description	Qty
23A00W669	24 Square Well, 10mL, V-Bottom, 41.9mm Height, Sterile	50
23A00W683	Elution Plate, 96 Square Well, 0.5mL, V-Bottom, 15.1mm Height, Sterile	50
23A00W693	96 Square Well, 1.0mL, V-Bottom, 24.4mm Height, Sterile	100
23A00W713	96 Square Well, 2.2mL, V-Bottom, 44.2mm Height, Sterile	50

KingFisher™ Compatible Deep Well Plates - Non-Sterile

Cat. No.	Description	Qty
23A00W668	24 Square Well, 10mL, V-Bottom, 41.9mm Height	50
23A00W682	Elution Plate, 96 Square Well, 0.5mL, V-Bottom, 15.1mm Height	50
23A00W692	96 Square Well, 1.0mL, V-Bottom, 24.4mm Height	100
23A00W712	96 Square Well, 2.2mL, V-Bottom, 44.2mm Height	50

*KingFisher™ is a trademark of Thermo Fisher Scientific

KingFisher™ Compatible Tip Combs

- Low-binding, medical grade polypropylene
- Compatible with KingFisher™ & Accuris IsoPure™
- 8, 24, and 96 tip configurations
- Available sterile or non-sterile

labForce® tip combs are produced from low-binding, medical grade polypropylene and are made specifically to fit KingFisher™ Flex, Accuris IsoPure™ and a variety of other automation platforms.

Applications include DNA extraction, quantification, viral RNA extraction, immunoprecipitations, etc.



KingFisher™ Compatible Tip Combs - *Sterile*

Cat. No.	Description	Qty
23A00W741	Tip Comb, 8-Tip, V-Bottom, PP, for KingFisher™, Sterile, 240 Bags of 2	480
23A00W743	Tip Comb, 24-Tip, V-Bottom, PP, for KingFisher™, Sterile, 25 Bags of 2	50
23A00W745	Tip Comb, 96-Tip, V-Bottom, PP, for KingFisher™, Sterile, 50 Bags of 2	100

KingFisher™ Compatible Tip Combs - *Non-Sterile*

Cat. No.	Description	Qty
23A00W740	Tip Comb, 8-Tip, V-Bottom, PP, for KingFisher™, 240 Bags of 2	480
23A00W742	Tip Comb, 24-Tip, V-Bottom, PP, for KingFisher™, 25 Bags of 2	50
23A00W744	Tip Comb, 96-Tip, V-Bottom, PP, for KingFisher™, 50 Bags of 2	100

*KingFisher™ is a trademark of Thermo Fisher Scientific