

## AMIGO2 Polyclonal Antibody

Catalog number: 14076-1-AP

Size: 49 µg/150 µl

Source: Rabbit

Isotype: IgG


Synonyms:

AMIGO2; ALI1, Alivin 1,

AMIGO 2, AMIGO2,

Amphoterin induced protein 2,

DEGA



Y79 cells were subjected to SDS PAGE followed by western blot with 14076-1-AP (AMIGO2 antibody) at dilution of 1:500

### Applications

Tested applications:	ELISA, WB
Species specificity:	Human, Mouse, Rat; other species not tested.
Calculated AMIGO2 MW:	57.9 kDa
Observed AMIGO2 MW:	58 kDa
Positive WB detected in	Y79 cells, K-562 cells
Recommended dilution:	WB: 1:500-1:5000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation

### Immunogen information

Immunogen:	Ag5200
GenBank accession number:	BC047595
Gene ID (NCBI):	347902
Full name:	Adhesion molecule with Ig-like domain 2

### Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3. Store at -20°C.