

# QUICK BLUE PROTEIN STAIN

\*for research & development purposes only

IBI Scientific's **Quick Blue Protein Stain** is a comprehensive rapid, 1-step coomassie type stain. The proprietary formulation, using colloidal coomassie, is ideal for **rapid protein staining** in polyacrylamide gels. This is a **non-hazardous** product that does NOT utilize organic solvents or phosphoric acid. A completely **environmentally friendly** product!

**MODEL**  
**NUMBER:**

IB01034 --- 1 Liter



## ADVANTAGES:

Detect proteins down to 5 ng

Does not contain organic solvents or phosphoric acid

Can be used for mass spectrometry