

Carboxypeptidase A5 Polyclonal Antibody

Catalog number: 13731-1-AP

Size: 20 µg/150 µl


Source: Rabbit

Isotype: IgG


Synonyms:

CPA5; Carboxypeptidase A5,

CPA5, FLJ35969



PC-3 cells were subjected to SDS PAGE followed by western blot with 13731-1-AP(CPA5 antibody) at dilution of 1:500



Immunofluorescent analysis of PC-3 cells using 13731-1-AP(CPA5 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)

Applications

Tested applications:	ELISA, WB, IF
Species specificity:	Human,Mouse,Rat; other species not tested.
Calculated Carboxypeptidase A5 MW:	436aa,49 kDa
Observed Carboxypeptidase A5 MW:	49 kDa
Positive WB detected in	PC-3 cells
Positive IF detected in	PC-3 cells
Recommended dilution:	WB: 1:500-1:5000 IF: 1:20-1:200

Application key: WB = Western blotting, IHC= Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation

Immunogen information

Immunogen:	Ag4755
GenBank accession number:	BC039362
Gene ID (NCBI):	93979
Full name:	Carboxypeptidase A5

Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3. Store at -20°C.