

DAP3 Polyclonal Antibody

Catalog number: 10276-1-AP

Size: 52 µg/150 µl

Source: Rabbit

Isotype: IgG


Synonyms:

DAP3; DAP-3;

DKFZp686G12159; MGC126058;

MGC126059; MRP-S29;

MRPS29; bMRP-10



HeLa cells were subjected to SDS PAGE followed by western blot with 10276-1-AP (DAP3 Antibody) at dilution of 1:1000

Applications

Tested applications:	ELISA, WB
Species specificity:	Human, mouse; other species not tested.
Positive WB detected in:	HeLa cells, human brain tissue
Calculated DAP3 MW:	46kd
Observed DAP3 MW:	40-45kd
Recommended dilution:	WB: 1:500-1:5000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation, FC = Flow cytometry

Immunogen information

Immunogen:	Ag0274
GenBank accession number:	BF220193
Gene ID (NCBI):	7818
Full name:	Death associated protein 3

Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3. Store at -20°C.