

|               |  |
|---------------|--|
| Product Name: | ABS ET 3   |
| Description:  | <p>ABS ET Series</p> <p>Includes 6 Canisters</p> <p>½ CC Straw (1 Level Bulk): 1122</p> <p>Static Evaporation Rate (Liters/day): 0.13</p> <p>LN2 Capacity: 3.6</p> <p>Normal Working Time: 17 Days</p> <p>Neck Diameter (mm): 55.4</p> <p>Exterior Height (mm): 406.4</p> <p>Exterior Diameter (mm): 220.9</p> <p>Canister Height (mm): 127</p> <p>Canister Diameter (mm): 41.9</p> <p>Designed for small capacity storage</p> <p>Durable aluminum tank with powder coated exterior</p> <p>High strength neck tube reduces liquid nitrogen loss</p> <p>Lid lock provision</p> <p>Low liquid nitrogen consumption</p> <p>Liquid nitrogen measuring rod for level validation</p> <p>One-year warranty on all parts and labor</p> <p>Three-year vacuum warranty</p> |