

HRH2 Polyclonal Antibody

Catalog number: 13414-1-AP


Size: 62 µg/150 µl

Source: Rabbit

Isotype: IgG

Synonyms:

HRH2; H2R



SH-SY5Y cells were subjected to SDS PAGE followed by western blot with 13414-1-AP (HRH2 Antibody) at dilution of 1:600

Applications

Tested applications:	ELISA, WB
Species specificity:	Human, mouse, rat; other species not tested.
Positive WB detected in:	SH-SY5Y cells, human skeletal muscle tissue, mouse stomach tissue
Calculated HRH2 MW:	397aa, 45kd
Observed HRH2 MW:	59kd, 69kd
Recommended dilution:	WB: 1:200-1:2000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation, FC = Flow cytometry

Immunogen information

Immunogen:	Ag4068
GenBank accession number:	BC054510
Gene ID (NCBI):	3274
Full name:	Histamine receptor H2

Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3. Store at -20°C.