



# BIOFIL®

Document Number

Revision

Page

Date

PMT001


A2

1 of 2

01-28-2021

## PRODUCT DATA SHEET

## ZEROTIP® Pipette Micro Tips



### ZEROTIP® Pipette Micro Tips

**Description:** Pipette tips, super low retention, natural, Dnase/Rnase-free, Non-pyrogenic, with Filter

**Purpose:** Use for liquid transfer

#### Materials

Tips: PP (Polypropylene)	Color: Natural
Rack: PP (Polypropylene)	Color: Color-coded
Filter: PO (Polyolefin)	Color: White


#### Features

- JET ZEROTIP® are manufactured from super clear high quality Polypropylene
- The surfaces of the tips are produced through a special process, this process makes the tip inner surface become super hydrophobic
- Preferred accessories for most brand micropipettor
- Tips with filter are also available
- Color-coded racks for easy volume identification
- Every rack or case is printed with lot No. for quality traceability
- DNase/RNase-free and non-pyrogenic
- Sterilized by irradiation SAL10<sup>-6</sup> (ISO11137)
- Shelf life: 3 years after month of production
- Manufactured in a Class 100,000 cleanroom environment
- Manufactured under ISO13485 and ISO9001 quality management system



STERILE | R

Dimensions ( Unit: mm)



### ZEROTIP® Pipette Micro Tips

Cat. No.	Volume	Tip Color	Surface Type	Filter	Sterilized	Packaging Spec.
PMT110010	10µL	Natural	Low retention	Yes	No	1000 tips/bag, 10 bags/case
PMT111010					Yes	
PMT550010					No	96 tips/rack, 10 racks/box, 2 boxes/case
PMT252010					Yes	
PMT711010					Yes	96 tips/bag, 20 bags/case
PMT130010	10µL, Long	Natural	Low retention	Yes	No	1000 tips/bag, 10 bags/case
PMT131010					Yes	
PMT232010					No	96 tips/rack, 10 racks/box, 2 boxes/case
PMT233010					Yes	
PMT731010					Yes	96 tips/bag, 20 bags/case

**BIOFIL<sup>®</sup>**

Document Number

Revision

Page

Date

PMT001

A2

2 of 2

01-28-2021

**PRODUCT DATA SHEET****ZEROTIP<sup>®</sup> Pipette Micro Tips****ZEROTIP<sup>®</sup> Pipette Micro Tips**

Cat. No.	Volume	Tip Color	Surface Type	Filter	Sterilized	Packaging Spec.
PMT110020	20µL	Natural	Low retention	Yes	No	1000 tips/bag, 10 bags/case
PMT111020					Yes	
PMT250020					No	96 tips/rack, 10 racks/box, 2 boxes/case
PMT252020					Yes	
PMT711020					Yes	
PMT110100	100µL	Natural	Low retention	Yes	No	1000 tips/bag, 10 bags/case
PMT111100					Yes	
PMT250100					No	96 tips/rack, 10 racks/box, 2 boxes/case
PMT252100					Yes	
PMT711100					Yes	
PMT012200	200µL	Natural	Low retention	Yes	No	1000 tips/bag, 10 bags/case
PMT111200					Yes	
PMT230200					No	96 tips/rack, 10 racks/box, 2 boxes/case
PMT231200					Yes	
PMT711200					Yes	
PMT130300	300µL	Natural	Low retention	Yes	No	1000 tips/bag, 10 bags/case
PMT131300					Yes	
PMT232300					No	96 tips/rack, 10 racks/box, 2 boxes/case
PMT233300					Yes	
PMT731300					Yes	
PMT110000	1000µL	Natural	Low retention	Yes	No	1000 tips/bag, 10 bags/case
PMT111000					Yes	
PMT550000					No	96 tips/rack, 10 racks/box, 2 boxes/case
PMT252000					Yes	
PMT711000					Yes	
PMT170000	1000µL, Long	Natural	Low retention	Yes	No	1000 tips/bag, 10 bags/case
PMT171000					Yes	
PMT370000					No	96 tips/rack, 10 racks/box, 2 boxes/case
PMT371000					Yes	

Drawn		Approved	
Date	28 Jan 2021	Date	28 Jan 2021
Name	Alex Fang	Name	James