

Hello “Tom”



The Thomas Optimized Mobile Cart from Thomas Scientific.
Specifically designed to bring rise to Personal Protective Equipment usage.

Thomas Optimized Mobile Cart (Tom)

Ergonomic design for optimal productivity

Upper module organizes critical inventory items at staff's fingertips for instant access with no bending, stooping or wasted reach. Gloves, masks, wipes, media, sterile packs, swabs, instruments, and more, all within easy reach.

Full-width, easy-grip, ergonomic drawer pulls. Dual-roller glides require minimal effort to open/close fully loaded drawers

Adjustable drawer dividers keep contents organized

Super-durable thermoplastic polymer construction ensures longevity:

- **Easily Sanitized** with common disinfectants and bleach solutions
- **Inert** - will not support germ/bacteria/fungal growth
- **UV-Stable** - colors will not fade, delaminate or deteriorate
- **Impact and Scratch-Resistant**
- **Waterproof** - safe to pressure-clean and store outdoors
- **FDA-Compliant** - Reg. 21, Chapter 1, part 177.2510

Load rated @ 400 lbs



Advanced design features

Adjustable dividers organize gloves, tissues, supplies, instruments, and critical workflow for easy access

Radiused 3mm safety edges to protect fingers, hands, walls and equipment

Ergonomic ADA pulls @34" height with soft grip for easy push-pull-turn cart navigation

Convenient prep/posting worksurface @ 34" high

Drawer locks are Keyed (standard), Digital Keypad or RFID access

Soft-close, non-drift drawers stay closed when cart is in motion

Safe & easy movement

Smooth, quiet casters with sealed precision bearings for easy rolling with minimal start-stop force

“Total lock” brake for effortless positioning and mobility of Mac



Website
www.ThomasSci.com



Email
sales@thomassci.com




Phone
800.345.2100

Follow Us:

Thomas Scientific, LLC Family of Companies



Specifications



Website
www.ThomasSci.com



Email
sales@thomassci.com



Phone
800.345.2100

Follow Us: [f](#) [@](#) [t](#) [in](#)



Thomas Scientific, LLC Family of Companies