

HOXC10 Polyclonal Antibody

Catalog number: 12025-1-AP

Size: 36 µg/150 µl

Source: Rabbit


Isotype: IgG

Synonyms:

HOXC10; HOXC10, homeobox
C10, HOX3I, Homeobox protein

Hox C10, Homeobox protein

Hox 3I



mouse heart tissue were subjected to SDS PAGE followed by western blot with 12025-1-AP(HOXC10 antibody) at dilution of 1:500

Applications

Tested applications:	ELISA, WB
Species specificity:	Human, mouse, rat; other species not tested.
Calculated HOXC10 MW:	342aa,38 kDa
Observed HOXC10 MW:	50kd
Positive WB detected in	Mouse heart tissue, A549 cells
Recommended dilution:	WB: 1:500-1:5000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation

Immunogen information

Immunogen:	Ag2655
GenBank accession number:	BC001293
Gene ID (NCBI):	3226
Full name:	Homeobox C10

Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3. Store at -20°C.