

PTS Polyclonal Antibody

Catalog number: 12150-1-AP

Size: 33 µg/150 µl


Source: Rabbit

Isotype: IgG


Synonyms:

PTS; FLJ97081, PTP synthase,

PTPS, PTS



Immunohistochemical of paraffin-embedded human gliomas using 12150-1-AP (PTS antibody) at dilution of 1:50 (under 10x lens)



human brain tissue were subjected to SDS PAGE followed by western blot with 12150-1-AP (PTS antibody) at dilution of 1:200

Applications

Tested applications:	ELISA, WB, IHC
Species specificity:	Human, Mouse, Rat; other species not tested.
Calculated PTS MW:	145aa, 16 kDa
Observed PTS MW:	16 kDa
Positive WB detected in:	Human brain tissue
Positive IHC detected in:	Human gliomas tissue
Recommended dilution:	WB: 1:200-1:2000 IHC: 1:20-1:200

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation

Immunogen information

Immunogen:	Ag2798
GenBank accession number:	BC009686
Gene ID (NCBI):	5805
Full name:	6-pyruvoyltetrahydropterin synthase

Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.1% sodium azide and 50% glycerol pH 7.3. Store at -20°C.