

ARMER Polyclonal Antibody

Catalog number: 12988-1-AP


Size: 42 µg/150 µl

Source: Rabbit

Isotype: IgG

Synonyms:

ARL6IP1; Aip 1, AIP1, ARL 6
interacting protein 1, ARL6IP,
ARL6IP1, ARMER, KIAA0069



mouse lung tissue were subjected to SDS PAGE followed by western blot with 12988-1-AP (ARL6IP1 antibody) at dilution of 1:300

Applications

Tested applications:	ELISA, WB
Species specificity:	Human, Mouse, Rat; other species not tested.
Calculated ARMER MW:	203aa, 23kd
Observed ARMER MW:	22-24kd
Positive WB detected in:	Mouse lung tissue, K-562 cells
Recommended dilution:	WB: 1:200-1:2000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation

Immunogen information

Immunogen:	Ag3639
GenBank accession number:	BC010281
Gene ID (NCBI):	23204
Full name:	ADP-ribosylation factor-like 6 interacting protein 1

Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3. Store at -20°C.