

VPS4A Polyclonal Antibody

Catalog number: 14272-1-AP

Size: 20 µg/150 µl

Source: Rabbit


Isotype: IgG

Synonyms:

VPS4A; FLJ22197, hVPS4,

Protein SKD2, SKD1, SKD2,

VPS4, VPS4 1, VPS4A



SMMC-7721 cells were subjected to SDS PAGE followed by western blot with 14272-1-AP(VPS4A Antibody) at dilution of 1:300

Applications

Tested applications:	ELISA, WB
Species specificity:	Human,Mouse,Rat; other species not tested.
Calculated VPS4A MW:	49 kDa
Observed VPS4A MW:	48 kDa
Positive WB detected in	SMMC-7721 cells
Recommended dilution:	WB: 1:200-1:1000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation

Immunogen information

Immunogen:	Ag5592
GenBank accession number:	BC047932
Gene ID (NCBI):	27183
Full name:	Vacuolar protein sorting 4 homolog A (S. cerevisiae)

Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.1% sodium azide and 50% glycerol pH 7.3. Store at -20°C.