

## Synaptotagmin-11 Polyclonal Antibody

Catalog number: 12031-1-AP

Size: 20 µg/150 µl

Source: Rabbit

Isotype: IgG

Synonyms:


SYT11; DKFZp781D015,

KIAA0080, Synaptotagmin 11,

synaptotagmin XI,

Synaptotagmin-11, SYT11,

SYT12, SytXI



mouse skeletal muscle tissue were subjected to SDS PAGE followed by western blot with 12031-1-AP(SYT11 antibody) at dilution of 1:300

### Applications

Tested applications: ELISA, WB

Species specificity: Human, Mouse, Rat; other species not tested.

Calculated Synaptotagmin-11 431aa, 48 kDa

MW:

Observed Synaptotagmin-11 96kd

MW:

Positive WB detected in Mouse skeletal muscle tissue

Recommended dilution: WB: 1:500-1:5000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation

### Immunogen information

Immunogen: Ag2663

GenBank accession number: BC013690

Gene ID (NCBI): 23208

Full name: Synaptotagmin XI

### Product information

Purification method: Antigen affinity purification

Storage: PBS with 0.02% sodium azide and 50% glycerol pH 7.3. Store at -20°C.