

FGFR1OP2 Polyclonal Antibody

Catalog number: 11605-1-AP

Size: 25 µg/150 µl

Source: Rabbit

Isotype: IgG


Synonyms:

FGFR1OP2; FGFR1OP2,

DKFZp564O1863, FGFR1

oncogene partner 2, HSPC123

like



mouse thymus tissue were subjected to SDS PAGE followed by western blot with 11605-1-AP(FGFR1OP2 antibody) at dilution of 1:500

Applications

Tested applications:	ELISA, WB
Species specificity:	Human, mouse, rat; other species not tested.
Calculated FGFR1OP2 MW:	29 kDa
Observed FGFR1OP2 MW:	25kd
Positive WB detected in	Mouse thymus tissue
Recommended dilution:	WB: 1:500-1:5000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation

Immunogen information

Immunogen:	Ag2153
GenBank accession number:	BC032143
Gene ID (NCBI):	26127
Full name:	FGFR1 oncogene partner 2

Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.1% sodium azide and 50% glycerol pH 7.3. Store at -20°C.