

HK3 Polyclonal Antibody

Catalog number: 13333-1-AP


Size: 28 µg/150 µl

Source: Rabbit


Isotype: IgG

Synonyms:


HK3; Hexokinase 3, hexokinase 3 (white cell), Hexokinase type III, HK III, HK3, HKIII, HXK3



Immunohistochemical of paraffin-embedded human lung cancer using 13333-1-AP(HK3 antibody) at dilution of 1:100 (under 10x lens)



Immunohistochemical of paraffin-embedded human lung cancer using 13333-1-AP(HK3 antibody) at dilution of 1:100 (under 40x lens)



Transfected HEK-293 cells were subjected to SDS PAGE followed by western blot with 13333-1-AP(HK3 Antibody) at dilution of 1:1000

Applications

Tested applications:	ELISA, IHC, WB
Species specificity:	Human, Mouse, Rat; other species not tested.
Calculated HK3 MW:	923aa, 99 kDa
Observed HK3 MW:	
Positive WB detected in	Transfected HEK-293 cells
Positive IHC detected in	Human lung cancer tissue
Recommended dilution:	WB: 1:500-1:5000 IHC: 1:20-1:200

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation

Immunogen information

Immunogen:	Ag4148
GenBank accession number:	BC028129
Gene ID (NCBI):	3101
Full name:	Hexokinase 3 (white cell)

Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.1% sodium azide and 50% glycerol pH 7.3. Store at -20°C.