

FBXO8 Polyclonal Antibody

Catalog number: 11903-1-AP


Size: 65 µg/150 µl

Source: Rabbit

Isotype: IgG

Synonyms:

FBXO8; DC10, F box only
protein 8, F box protein 8, F
box/SEC7 protein FBS, FBS,
FBX8, FBXO8,
UNQ1877/PRO4320



human placenta tissue were
subjected to SDS PAGE followed by
western blot with 11903-1-AP(FBXO8
antibody) at dilution of 1:500

Applications

Tested applications:	ELISA, WB
Species specificity:	Human, Mouse, Rat; other species not tested.
Calculated FBXO8 MW:	37 kDa or 32 kDa
Observed FBXO8 MW:	45-50 kDa
Positive WB detected in	Human placenta tissue, HeLa cells, mouse testis tissue
Recommended dilution:	WB: 1:500-1:5000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation

Immunogen information

Immunogen:	Ag2539
GenBank accession number:	BC040456
Gene ID (NCBI):	26269
Full name:	F-box protein 8

Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.1% sodium azide and 50% glycerol pH 7.3. Store at -20°C.