

Sharp™ Pipette Starter Kit

Sharp™ pipettes are now available in a convenient starter kit for a great price.

- Fully autoclavable
- Reduced plunger forces
- The leader in accuracy

Each starter kit includes:

- (4) Ultra EZpette™ pipettes (0.5–10µL, 2–20µL, 20–200µL and 100–1000µL)
- (1) six position carousel stand
- (3) full racks each of pipette tips (10µl, 200µl, 1250µl)



Thomas No.	Description	Quantity
20A00C937	Sharp™ Pipette Starter Kit	Each

MAY WE ALSO SUGGEST:

SHARP™ Single and Multi-Channel Pipettors

SHARP™ Pipettors provide both accurate and precise liquid dispensing in a comfortable ergonomic design, dual volume setting and a volume lock to avoid volume changes during routine pipetting. The ejector “klik-to-fix” system easily adjusts tip ejection.

SHARP™ Pipettors are fully autoclavable. Eight models cover the volume range from 0.1µL to 10,000µL for single channel, and 1µL to 300µL for multi-channel. A calibration tool is included with each unit for easy in-lab recalibration. Pipettes are individually tested and supplied with a certificate of quality.



For product selection, visit: ThomasSci.com

For Single channel, search: [SHARP™ Single Channel](#)

For Multi-Channel, search: [SHARP™ Multi Channel](#)