

Repeat Pipetting is *EASIER* with

ali-Q 2TM aliquoting pipet controller


Variable aliquot
Speed (VS) model
for all your
application needs




• aliquot speed dial



• aliquot button



• aliquot volume dial



The ali-Q 2 Series:

Variable aliquoting speed model

- ali-Q 2 VS

VistaLab
Catalog #

2100-1005

max. aliquot
volume

5 mL

Fixed aliquoting speed models

- ali-Q 2 (high speed)
- ali-Q 2 LS (Low Speed)

2100-0005

2100-0003

5 mL

3 mL

- **Repeat pipetting** / multidispensing with less fatigue
- Highly accurate and reproducible aliquots with one touch
- Use Serological Pipets but no need to watch the lines
- In-lab calibration with on-boarded software

Wobble-not™

Low Insertion Force Serological Pipet

**Unique, patented
No Drip Design!**



Do any of your other serological pipets look like this?

2-tier step design

- = 2 zones of contact in your pipet controller = Less force to insert
- = Less "wobble" = Max stability
- = Better seal = Less dripping = Better accuracy

LOSE THE WOBBLE

**LOSE
THE
DRIP**



Ordering Information

| Pipet Volume | VistaLab Catalog # | Quantity per case |
|--------------|--------------------|-------------------|
|--------------|--------------------|-------------------|

Standard length

Individually wrapped (paper/plastic), Sterile

| | | |
|-------|-----------|-----|
| 5mL | 4090-0005 | 200 |
| 10mL | 4090-0010 | 200 |
| 25mL | 4090-0025 | 200 |
| 50mL | 4090-0050 | 100 |
| 100mL | 4090-0100 | 50 |

Individually wrapped (plastic/plastic), Sterile

| | | |
|-------|-----------|-----|
| 5mL | 4190-0005 | 200 |
| 10mL | 4190-0010 | 200 |
| 25mL | 4190-0025 | 200 |
| 50mL | 4190-0050 | 100 |
| 100mL | 4190-0100 | 50 |

Bulk packaged (resealable plastic bag), Sterile



| | | |
|-------|-----------|-----|
| 5mL | 4290-0005 | 500 |
| 10mL | 4290-0010 | 400 |
| 25mL | 4290-0025 | 150 |
| 50mL | 4290-0050 | 100 |
| 100mL | 4290-0100 | 60 |

Shorty length

Individually wrapped (paper/plastic), Sterile

| | | |
|------|-----------|-----|
| 10mL | 5090-0010 | 200 |
| 25mL | 5090-0025 | 100 |

Why it works



VistaLab™
TECHNOLOGIES

Think Differently. Feel Better.