

GLTP Polyclonal Antibody

Catalog number: 10850-1-AP


Size: 37 µg/150 µl

Source: Rabbit

Isotype: IgG

Synonyms:

GLTP; GLTP, glycolipid transfer protein



mouse brain tissue were subjected to SDS PAGE followed by western blot with 10850-1-AP (GLTP antibody) at dilution of 1:500

Applications

Tested applications:	ELISA, WB
Species specificity:	Human, Mouse, Rat; other species not tested.
Calculated GLTP MW:	24 kDa
Observed GLTP MW:	24 kDa
Positive WB detected in	Mouse brain tissue, A2780 cells
Recommended dilution:	WB: 1:200-1:2000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation

Immunogen information

Immunogen:	Ag1293
GenBank accession number:	BC009932
Gene ID (NCBI):	51228
Full name:	Glycolipid transfer protein

Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.1% sodium azide and 50% glycerol pH 7.3. Store at -20°C.