

TFDP1 Polyclonal Antibody

Catalog number: 11043-1-AP

Size: 37 µg/150 µl

Source: Rabbit

Isotype: IgG


Synonyms:

TFDP1; Dp 1, DP1, DRTF1,

DRTF1 polypeptide 1, E2F

dimerization partner 1, TFDP1,

transcription factor Dp 1



Immunohistochemical of paraffin-embedded human lymphoma using 11043-1-AP (TFDP1 antibody) at dilution of 1:50 (under 10x lens)

Applications

Tested applications:	ELISA, IHC
Species specificity:	Human, Mouse, Rat; other species not tested.
Calculated TFDP1 MW:	46kd
Observed TFDP1 MW:	
Positive IHC detected in	Human lymphoma tissue
Recommended dilution:	IHC: 1:25 - 1:100

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation

Immunogen information

Immunogen:	Ag1511
GenBank accession number:	BC011685
Gene ID (NCBI):	7027
Full name:	Transcription factor Dp-1

Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.1% sodium azide and 50% glycerol pH 7.3. Store at -20°C.