

For Research Use Only

IFT80 POLYCLONAL ANTIBODY



Catalog Number: 25230-1-AP

Basic Information

Catalog Number:
25230-1-AP

Size:
66 µg/150 µl

Source:
Rabbit

Isotype:
IgG

Purification Method:
Antigen affinity purification

Immunogen Catalog Number:
AG18650

GenBank Accession Number:
BC101494

GeneID (NCBI):
57560

Full Name:
intraflagellar transport 80 homolog
(Chlamydomonas)

Calculated MW:
777aa, 88 kDa

Observed MW:
88-90 kDa

Recommended Dilutions:
WB 1:500-1:5000

Applications

Tested Applications:
ELISA, WB

Species Specificity:
Human, Mouse

Positive Controls:

WB : NIH/3T3 cells; mouse brain tissue, mouse kidney tissue

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
--------	-----------	---------	-------------

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.1% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



NIH/3T3 cells were subjected to SDS PAGE followed by western blot with 25230-1-AP (IFT80 Antibody) at dilution of 1:1000