

CAMK1D Polyclonal Antibody

Catalog number: 13613-1-AP

Size: 20 µg/150 µl


Source: Rabbit

Isotype: IgG

Synonyms:

CAMK1D; CKLiK; CaM-K1;

CaMKID



mouse thymus tissue were subjected to SDS PAGE followed by western blot with 13613-1-AP(CAMK1D Antibody) at dilution of 1:2000

Applications

Tested applications:	ELISA, WB
Species specificity:	Human, mouse, rat; other species not tested.
Positive WB detected in:	Mouse thymus tissue
Calculated CAMK1D MW:	357aa,40kd
Observed CAMK1D MW:	40kd
Recommended dilution:	WB: 1:1000-1:10000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation, FC = Flow cytometry

Immunogen information

Immunogen:	Ag4517
GenBank accession number:	BC035745
Gene ID (NCBI):	57118
Full name:	Calcium/calmodulin-dependent protein kinase ID

Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3. Store at -20°C.