

## TRIM22 Polyclonal Antibody

Catalog number: 13744-1-AP


Size: 45 µg/150 µl

Source: Rabbit

Isotype: IgG

Synonyms:

TRIM22; GPSTAF50, RING  
finger protein 94, RNF94, Staf  
50, STAF50, TRIM22, tripartite  
motif containing 22



mouse thymus tissue were subjected to SDS PAGE followed by western blot with 13744-1-AP( TRIM22 Antibody) at dilution of 1:300

### Applications

Tested applications:	ELISA, WB
Species specificity:	Human, Mouse; other species not tested.
Calculated TRIM22 MW:	498aa,57 kDa
Observed TRIM22 MW:	50 kDa
Positive WB detected in	Mouse thymus tissue
Recommended dilution:	WB: 1:200-1:1000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation

### Immunogen information

Immunogen:	Ag4683
GenBank accession number:	BC035582
Gene ID (NCBI):	10346
Full name:	Tripartite motif-containing 22

### Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.1% sodium azide and 50% glycerol pH 7.3. Store at -20°C.