

AAIMI Level 2 Isolation Gowns

AAIMI Levels 2 Isolation Gowns provide barrier protection to help reduce the risk of contamination and transmission of infectious organisms.

Made of a triple-layer SMS fabric, they are ideal for low-risk procedures, such as drawing blood from veins and suturing small wounds.

Features:

- Size: 30" x 46"
- Yellow SMS material that meets AAMI PB70 Level 2 testing
- Strong, welded seams
- Waist belt and elastic cuffs
- Full back with velcro tab neck tie
- Not made with natural rubber latex
- One size fits all



Thomas No.	Description	Quantity
20A00H355	Isolation Gown, AAIMI Level 2, Yellow SMS, Size X-Large	100/Case



Disclaimer: Pricing and offers subject to change without notice.