

## DOK6 Polyclonal Antibody

Catalog number: 14108-1-AP

Size: 30 µg/150 µl


Source: Rabbit

Isotype: IgG

Synonyms:

DOK6; DOK5L; HsT3226;

MGC20785



mouse brain tissue were subjected to SDS PAGE followed by western blot with 14108-1-AP(DOK6 Antibody) at dilution of 1:1000

### Applications

Tested applications:	ELISA, WB
Species specificity:	Human, mouse, rat; other species not tested.
Positive WB detected in:	Mouse brain tissue
Calculated DOK6 MW:	38kd
Observed DOK6 MW:	38kd
Recommended dilution:	WB: 1:500-1:5000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation, FC = Flow cytometry

### Immunogen information

Immunogen:	Ag5254
GenBank accession number:	BC008583
Gene ID (NCBI):	220164
Full name:	Docking protein 6

### Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3. Store at -20°C.