

UBL3 Polyclonal Antibody

Catalog number: 14100-1-AP

Size: 43 µg/150 µl


Source: Rabbit

Isotype: IgG

Synonyms:

UBL3; DKFZp434K151;

FLJ32018; HCG-1; PNSC1



human brain tissue were subjected to SDS PAGE followed by western blot with 14100-1-AP(UBL3 Antibody) at dilution of 1:1500

Applications

Tested applications:	ELISA, WB
Species specificity:	Human, mouse, rat; other species not tested.
Positive WB detected in:	Human brain tissue
Calculated UBL3 MW:	117aa,13kd
Observed UBL3 MW:	13kd
Recommended dilution:	WB: 1:500-1:5000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation, FC = Flow cytometry

Immunogen information

Immunogen:	Ag5233
GenBank accession number:	BC059385
Gene ID (NCBI):	5412
Full name:	Ubiquitin-like 3

Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3. Store at -20°C.