

STAC Polyclonal Antibody

Catalog number: 11480-1-AP


Size: 30 µg/150 µl

Source: Rabbit

Isotype: IgG

Synonyms:

STAC; FLJ32331; STAC1



PC-3 cells were subjected to SDS PAGE followed by western blot with 11480-1-AP (STAC Antibody) at dilution of 1:2000

Applications

Tested applications:	ELISA, WB
Species specificity:	Human, mouse, rat; other species not tested.
Positive WB detected in:	PC-3 cells, human liver tissue
Calculated STAC MW:	402aa, 45kd
Observed STAC MW:	45kd
Recommended dilution:	WB: 1:1000-1:10000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP =

Immunoprecipitation, FC = Flow cytometry

Immunogen information

Immunogen:	Ag2059
GenBank accession number:	BC020221
Gene ID (NCBI):	6769
Full name:	SH3 and cysteine rich domain

Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.1% sodium azide and 50% glycerol pH 7.3. Store at -20°C.