

## B4GALT3 Polyclonal Antibody

Catalog number: 11041-1-AP


Size: 30 µg/150 µl

Source: Rabbit

Isotype: IgG

Synonyms:

B4GALT3; beta4Gal-T3



Immunohistochemistry of paraffin-embedded human breast cancer using 11041-1-AP (B4GALT3 Antibody) at Dilution 1:50 (under 10x lens)

### Applications

Tested applications:	ELISA, IHC
Species specificity:	Human, mouse, rat; other species not tested.
Calculated B4GALT3 MW:	44kd
Observed B4GALT3 MW:	
Positive IHC detected in:	Human breast cancer
Recommended dilution:	IHC: 1:20-1:200

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation, FC = Flow cytometry

### Immunogen information

Immunogen:	Ag1487
GenBank accession number:	BC009985
Gene ID (NCBI):	8703
Full name:	UDP-Gal:betaGlcNAc beta 1,4- galactosyltransferase, polypeptide 3

### Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.1% sodium azide and 50% glycerol pH 7.3. Store at -20°C.