

TRA Polyclonal Antibody

Catalog number: 11890-1-AP


Size: 49 µg/150 µl

Source: Rabbit


Isotype: IgG

Synonyms:

TRA-; FLJ22602, T cell receptor
alpha locus, TCRA, TCRD, TRA,
TRA-



Immunohistochemistry of paraffin-embedded human tonsillitis tissue slide using 11890-1-AP (TRA-Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded human tonsillitis tissue slide using 11890-1-AP (TRA-Antibody) at dilution of 1:50 (under 40x lens)

Applications

Tested applications:	ELISA, IHC
Species specificity:	Human; other species not tested.
Calculated TRA MW:	202aa,23 kDa
Observed TRA MW:	
Positive IHC detected in	Human tonsillitis tissue
Recommended dilution:	IHC: 1:20-1:200

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation

Immunogen information

Immunogen:	Ag2498
GenBank accession number:	BC020840
Gene ID (NCBI):	6955
Full name:	T cell receptor alpha locus

Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.1% sodium azide and 50% glycerol pH 7.3. Store at -20°C.