

Captair 321

Ventilated balance safety enclosure (VBSE)

The Captair Smart 321 protects users from both chemical vapors and airborne particulates. Configured with HEPA and/or carbon filtration, it's built to handle powder weighing applications with precision, stable enough to weigh to six decimal places without balance interference. No ductwork, no HVAC connection, and Smart Technology that alerts you the moment anything is off. Safety validated to your exact chemistry through Erlab's ValiQuest** process before it ever ships.



- All units will sit on existing bench in a lab
- Available in 14 business days

Powder & Liquid Weighing Station - HEPA & CARBON**

**VQ required as liquid chemicals will be used.

Order Part Number	Includes	QTY
26A00B554	321 Filter Housing	1
	321 Enclosure	1
	321 Resin Worksurface	1
	Hepa Filter	1
	Molecode- S	1
	Carbon Filter	1

New Bundled Solutions:

General Powder Weighing Station - HEPA ONLY*

*NO VQ is required, sales rep must confirm no liquid usage is occurring.

Order Part Number	Includes	QTY
26A00B552	321 Filter Housing	1
	321 Enclosure	1
	321 Resin Worksurface	1
	Hepa Filter	1
	Pre-Filter M	1

Carbon Only (BE+)**

**VQ required as liquid chemicals will be used.

Order Part Number	Includes	QTY
26A00B553	321 Filter Housing	1
	321 Enclosure	1
	321 Resin Worksurface	1
	Molecode- S	1
	Carbon Filters	2

Bundle offer expires 12/31/2026

