

# Sterile Swabs

steripackgroup.com





**3" Sterile Polyester Nasal Swabs used for Sample Collection to help combat Covid-19.**

# Sterile Swabs

Each swab is individually packaged and labeled in SteriPack's state-of-the-art ISO Class 8 Cleanroom manufacturing facilities located across North America, Europe & Asia. The shorter swab allows for a much more comfortable user experience when compared to that of a long Nasopharyngeal (NP) swab. Once the sample has been taken the swab can be placed straight in a tube rather than broken and then placed. Studies have shown it can be transported in a dry tube, removing the need for expensive Viral transport media (VTM).

## Specification SteriPack 3" Polyester Nasal Swab

**New** stronger binder  
for pool testing

	60564 	60577 
Sterile packed	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
100% polyester bud	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Polypropylene stick	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Gamma sterilized	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
CE Marked	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Bud weight (10 swabs)	0.31 grams	0.20 grams
Bud diameter	0.24 in. (6.1 mm)	0.2 in. (5.08 mm)
Swab length	2.87 in. (73 mm)	2.835 in. (72 mm)
Stick diameter	0.098 in. (2.5 mm)	0.1024 in. (2.6 mm)
Stick length	2.756 in. (70 mm)	2.756 in. (70 mm)
Bud fibers keep integrity and adhesion to stick under normal use		

SteriPack USA Ltd LLC  
4255 South Pipkin Road,  
Lakeland, Florida, 33811, USA

**Order Now**  
**+1 863 648 2333**  
swabs@steripackgroup.com  
steripackgroup.com